

SARS-CoV-2 Diagnostic Solutions; inqaba biotec



inqaba biotec

Africa's Genomics Company




Dr Aron Abera (Technical Support Manager)

aron.abera@inqababiotec.co.za

www.Inqababiotec.co.za

About inqaba




SARS-CoV-2


Severe Acute Respiratory Syndrome CoronaVirus 2

- RNA Virus
- Cause of COVID-19 disease

Relative positions of amplicon targets on the SARS coronavirus and the 2019 novel coronavirus genome



Corman et al 2020




Approved Workflow of SARS-CoV2 Molecular Diagnosis




Collect Specimens



Extract RNA



Analysis



Specimen Collection

Challenges

- Safety regarding specimen collection, handling & transportation



Solution

- DNA/RNA Shield™ Collection Devices from Zymo
- Collect specimens using DNA/RNA Shield™ Collection Devices
- SARS-CoV-2 RNA is stable in DNA/RNA Shield™ for **28 days** at room temperature
- Samples in DNA/RNA Shield™ can be shipped at ambient temperature – no cold chain needed

RNA Extraction

Challenges

- Low throughput RNA extraction and RNA extraction not optimized for biological fluids/ clinical samples

Solution



- RNA Extraction and Purification Quick-DNA/RNA Viral MagBead Kit from Zymo Research
 - Automated and Manual RNA extraction
 - Optimized for Biological fluids (Upper and lower respiratory specimen)
 - Compatible with UTM/VTM, Saline, PBS and liquid amies
- Monarch Total RNA Miniprep kit

RNA Extraction (Automation)



Primers and Probes

- Custom primer and probes for SARS-CoV-2 qPCR detection synthesized locally by inqaba biotec
- Primers and probes synthesized with speed, accuracy and purity
- Compatible probes synthesized for multiplexing several targets of SARS-CoV-2 and human internal control (e.g. FAM, HEX and others)



Facilitating Detection of SARS-CoV-2 Individual Reagents

RNA detection Reagents (NEB)

Challenges

- The inclusion of an up-front RNA purification step, adds to protocol time and cost per sample, and overall can reduce testing throughput

Solution



- Luna Universal Probe One-Step RT-qPCR kit
 - Hydrolysis probe-based detection of RNA
- LunaScript RT SuperMix kit
 - First strand synthesis in two-step RT-qPCR workflows
- WarmStart Colorimetric LAMP Master mix
 - For rapid reverse transcription and isothermal amplification of RNA targets, with a simple visible detection

Sequencing and COVID-19 (NEB)

- Illumina Sequencing (Reagents and Services)
 - Using the ARTIC SARS-CoV-2 primers (400bp) sample pooling and sequencing on the Illumina
 - Q5® Hot Start High-Fidelity DNA Polymerase
 - NEBNext® Ultra™ II End Repair / dA-Tailing Module
 - NEBNext Ultra II Ligation Module
 - NEBNext Quick Ligation Module
- PacBio HiFi Sequencing (Reagents and Services)
 - Long amplicon virus sequencing
- Nanopore Sequencing
 - Amplicon sequencing protocol for MinION for nCoV-2019,
 - Q5® Hot Start High-Fidelity DNA Polymerase
 - NEBNext® Ultra™ II End Repair / dA-Tailing Module
 - NEBNext Ultra II Ligation Module
 - NEBNext Quick Ligation Module

Allplex SARS-CoV-2 Assay (Seegene)

- Seegene, a company with wide range of solutions applicable in diverse markets from clinical diagnostics to basic science research




- Auto extraction System
 - Seegene STARlet/Nimbus
 - SEEPREP32
 - NucliSENS eayMAG
- Real Time PCR
 - CFX96 SW IVD v1.6 and Dx System v3.1
 - AB 7500 SW v2.05



Allplex SARS-CoV-2 Assay (Seegene)

- Four target genes assay
- Testing principle Real time PCR
- Multiplex Real-time PCR assay to detect 4 targets (RdRP, S and N gene specific for SARS CoV-2 and E gene for all of Sarbecovirus including SARS-CoV-2)
- Exogenous gene is used as an internal control (IC) to monitor the whole process.



Allplex SARS-CoV-2 Assay (Seegene)

- Sample Types
 - Nasopharyngeal/throat swab, Nasopharyngeal aspirate bronchoalveolar lavage and sputum
- Analytical Sensitivity
 - 50 copies/rxn for all targets
- Cut-offs
 - All analytes: ≤ 40 cycle
 - IC: ≤ 40 cycle

Allplex SARS-CoV-2 Assay (Seegene)

- Supplies and Accessories needed
 - RNA stabilization media
 - RNA extraction kit if not using NIMBUS automated extraction kit
 - Laboratory Consumables (Gloves, pipettes, plates, tips ...)



- Contents
 - Oligo Mix
 - PREMIX
 - Control (SARS2 PC)
 - Control IC
 - Rnase-free water

Test kits (50 and 100 reactions) and should be stored at -20° C


Allplex SARS-CoV-2 Assay (Seegene)



Interpretation of Results

	Case 1	Case 2	Case 3	Case 4	Case 5	Case 6	Case 7	Case 8	Case 9
IC (HEX)	+/-	+/-	+/-	+/-	+/-	+/-	+/-	+	-
E gene (FAM)	+	+	-	+	-	-	+	-	-
RdRP gene / S gene (Cal Red 610)	+	-	+	+	+	-	-	-	-
N gene (Quasar 670)	+	+	+	-	-	+	-	-	-
Result Interpretation	SARS-CoV-2 Detected						Presumptive positive	Negative	Invalid

Allplex SARS-CoV-2 Assay (Seegene)



APPLY RESULT

Well Info Positive Find Vertical Horizontal

Sample No	Patient Id	Well	Name	Type	FAM		Cal Red 610		Quasar 670		HEX		Auto Interpretation	Comment
					E gene	C(t)	RdRP/S gene	C(t)	N gene	C(t)	IC	C(t)		
<input type="checkbox"/>		E05		SAMPLE	+	36,07	+	36,37	+	35,89	+	22,38	SARS-CoV-2	
<input type="checkbox"/>		F05		SAMPLE	+	36,15	+	36,57	+	35,39	+	22,51	SARS-CoV-2	
<input type="checkbox"/>		G05		SAMPLE	+	24,33	+	24,14	+	25,11	+	23,24	SARS-CoV-2	
<input type="checkbox"/>		H05		SAMPLE	-	N/A	-	N/A	-	N/A	+	22,94	-	
<input type="checkbox"/>		A06		SAMPLE	-	N/A	-	N/A	+	36,95	+	22,72	SARS-CoV-2	
<input type="checkbox"/>		B06		SAMPLE	+	36,05	+	35,56	+	35,91	+	22,57	SARS-CoV-2	
<input type="checkbox"/>		C06		SAMPLE	-	N/A	+	34,51	+	36,46	+	23,18	SARS-CoV-2	
<input type="checkbox"/>		D06		SAMPLE	-	N/A	-	N/A	-	N/A	+	21,11	-	
<input type="checkbox"/>		E06		SAMPLE	+	35,62	+	37,16	+	36,49	+	22,92	SARS-CoV-2	
<input type="checkbox"/>		F06		SAMPLE	-	N/A	-	N/A	-	N/A	+	22,37	-	

THANK YOU.

