


# Roche SARS-CoV-2 diagnostic testing



*ASLM ECHO session #18, July 10<sup>th</sup> 2020*

---

*Dr. Wim van der Helm, MD - Healthcare Development EMEA-LATAM – Global Access Program*



*All statements made in this  
document are based on the current  
state of scientific literature*


July 2020

# Value of Diagnostic Tests in the New World of COVID-19 Pandemic




## The Problem

*In vitro diagnostic (IVD) tests are an essential service in the delivery of healthcare*



**Prior to pandemic IVD test results influence 70% of clinical decisions**




# Value of Diagnostic Tests in the New World of COVID-19 Pandemic



## Roche solution

*In vitro diagnostic (IVD) testing plays a vital role in testing who may have contracted COVID-19*



# What are we trying to find out?


## Testing Objectives

*“Should I advise isolation?”*

*“Could they still transmit the virus?”*


*“Have they previously been infected with the virus?”*

Is a person **currently infected** with SARS-CoV-2?




Has a person **previously been infected** with SARS-CoV-2?\*

*or*



# Human coronaviruses are zoonotic

*Originate in animals and can infect humans*



**Widest variety seen in bats**

Severe respiratory syndrome viruses

**7 known human coronaviruses**

Common cold viruses


**Mammals may be intermediate hosts**

Diagram Adapted from 1) and 3)

1. Zhu, N et al. (2020). *N Engl J Med* 2020;382:727-33; 2. Cui J, Li F, Shi ZL (2019). *Nature Reviews. Microbiology*. 17 (3): 181-92; 3. Su, S et al. (2016). *Trends in Microbiology*. 24 (6): 490-502; 4. Huynh J et al. (2012). *Journal of Virology*. 86 (23): 12816-25; 5. Masters PS (2006). *Advances in Virus Research*. Academic Press. 66: 193-292; 6. Vijgen, L. (2005). *Journal of Virology*. 79 (3): 1595-1604

# Human coronaviruses are zoonotic

*Originate in animals and can infect humans*



**Two previous outbreaks:**

- Severe acute respiratory syndrome – SARS-CoV (2002; China)
- Middle East respiratory syndrome – MERS-CoV (2012; Saudi Arabia)

Diagram Adapted from 1) and 3)

1. Zhu, N et al. (2020). *N Engl J Med* 2020;382:727-33; 2. Cui J, Li F, Shi ZL (2019). *Nature Reviews. Microbiology*. 17 (3): 181-92; 3. Su, S et al. (2016). *Trends in Microbiology*. 24 (6): 490-502; 4. Huynh J et al. (2012). *Journal of Virology*. 86 (23): 12816-25; 5. Masters PS (2006). *Advances in Virus Research*. Academic Press. 66: 193-292; 6. Vijgen, L. (2005). *Journal of Virology*. 79 (3): 1595-1604

# SARS-CoV-2 infection causes COVID-19

## Signs and Symptoms – an Overview<sup>1-7</sup>



### Transmission

Person-to-person via respiratory secretions; indirectly through contaminated surfaces



### Incubation

Range of 2–14 days (median – 5 days)



### Main Symptoms

Fever, respiratory symptoms, Anosmia, diarrhea, vomiting, headache, myalgia



### Clinical presentation

Asymptomatic infection, mild illness, pneumonia, or fatal disease



### Clinical progression

Can cause severe respiratory disease, especially in 65+ and multi-morbid patients



### Case fatality rate\*

Reported mortality rates vary from 0.8–15.5% (average 6.5%)

\* CFR is unreliable during an outbreak

(1) Chen N et al. (2020). Lancet. 2020;395:507–13; (2) Holshue ML et al. (2020). N Engl J Med. 382:929–36; (3) Huang C et al. (2020). Lancet 395:497–506; (4) Wang D et al. (2020). JAMA 323:1061–1069; (5) CDC. <https://www.cdc.gov/coronavirus/2019-nCoV/hcp/clinical-criteria.html> (6) Wu JT et al (2020). Lancet 395:689–97; (7) <https://informationisbeautiful.net/visualizations/covid-19-coronavirus-infographic-datapack/>

# COVID-19: risk factors for severe disease and death

*Accumulating evidence<sup>1</sup>*

Evidence is rapidly accumulating about risk factors for COVID-19 progression, with data suggesting that the following factors may be associated with increased risk of severe disease and/or mortality



**Older age**



**Male sex**



**Comorbidities**

(hypertension, diabetes, cardiovascular disease, chronic respiratory disease, cancer)



**Obesity**



**People who smoke**

The **relative importance of different underlying health conditions is unclear**, owing to:

Inadequate adjustment for important confounding factors, insufficient follow-up, and likely under-reporting of pre-existing conditions

Studies are mainly among those patients at the highest risk admitted to hospital with full testing and **might not apply to the general population**

# Coronaviruses

## *Virion morphology and structural proteins*

### Large enveloped RNA viruses (80–120 nm)

#### Lipid bilayer



Membrane protein (M) —

Envelope protein (E) —

Spike protein (S) —

#### Nucleocapsid

Multiple copies of the nucleocapsid (N protein) bound to the RNA genome




Masters PS (2006). Advances in Virus Research. Academic Press. 66: 193–292

Su, S et al. (2016). Trends in Microbiology. 24 (6): 490–502.

Paules CI et al. (2020). JAMA. 2020;323(8):707–708

# Following the infection path...

## *Viral transmission*



Virus enters  
the body via droplets



### Incubation

Range of 2–14 days  
(median – 5 days)




### Transmission

Person-to-person via respiratory  
secretions; indirectly through  
contaminated surfaces

# Following the infection path...

## *Viral reproduction in host cells of lung*


### Broad range of symptoms and clinical presentation



1. Li Q, et al. N Engl J Med (2020);doi.10.1056/NEJMoa2001316. 2. Backer JA, et al. Euro Surveill (2020);252000062.; 3. The World Health Organization <https://www.who.int/news-room/q-a-detail/q-a-coronaviruses>. Accessed 28 April 2020. 4. Zhou F, et al. Lancet (2020);395:1054–62. 5. Wu Y, et al. Lancet (2020);S2468–1253:30083. 1. Chen N, et al. Lancet (2020);395:507–513. PMID: 32007143; 2. Holshue ML, et al. N Engl J Med (2020);382:929–936. PMID: 32004427; 3. Huang C, et al. Lancet (2020);395:497–506. PMID: 31986264; 4. Wang D, et al. JAMA. (2020);323:1061–1069. PMID: 32031570; 5. Centers for Disease Control and Prevention. <https://www.cdc.gov/coronavirus/2019-nCoV/hcp/clinical-criteria.html>. Accessed March 4, 2020; 6. Wu JT, et al. Lancet (2020);395:689–697

# Following the infection path...


## *Stages of transient viral infections*



1. Li Q, et al. N Engl J Med (2020);doi.10.1056/NEJMoa2001316. 2. Backer JA, et al. Euro Surveill (2020);252000062.; 3. The World Health Organization <https://www.who.int/news-room/q-a-detail/q-a-coronaviruses>. Accessed 28 April 2020.  
4. Zhou F, et al. Lancet (2020);395:1054–62. 5. Wu Y, et al. Lancet (2020);S2468–1253:30083.

# Clinical stages of COVID-19

## Potential therapeutic approaches



# COVID-19 testing recommendations



## Criteria to guide evaluation and laboratory testing for COVID-19

<https://www.cdc.gov/coronavirus/2019-nCoV/hcp/clinical-criteria.html>



## Laboratory support for COVID-19 in the EU/EEA

<https://www.ecdc.europa.eu/en/novel-coronavirus/laboratory-support>



## Laboratory testing strategy recommendations for COVID-19

<https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/laboratory-guidance>


### In summary:

Nucleic Acid Amplification Test (NAAT) for diagnosis

Testing is needed to inform clinical decision making and public health policy

Prioritization of testing is needed when the availability of testing is limited

# Course of molecular and serological biomarkers\* SARS-CoV-2 infection




\*Illustrative only as seroconversion time not known as yet. Larger studies are required

Xiang F, et al. Clin Infect Dis (2020);ciaa461. doi: 10.1093/cid/ciaa461; 2. Zhao J, et al. Clin Infect Dis (2020);pii:ciaa344. doi:10.1093/cid/ciaa344; 3. Wölfel R, et al. Nature (2020);[epub ahead of print]. doi:10.1038/s41586-020-291-x. 4. Jin Y, et al. Int J Infect Dis (2020). 5. Liu W, et al. J Clin Microbiol (2020);pii:JCM.00461-20.doi:10.1128/JCM.00461-20. 6. Guo L, et al. Clin Infect Dis (2020); doi:10.1093/cid/ciaa310. 7. Zhang W, et al. Emerg Microb Infect;9:386-389.

# *Assessing the response to infection*

# Opportunities for detection of SARS-CoV-2 infection & recovery phases\*




\* Illustrative only as sero conversion time not known as yet. Larger studies are required

1. Xiang F, et al. Clin Infect Dis (2020);ciaa461. doi: 10.1093/cid/ciaa461; 2. Zhao J, et al. Clin Infect Dis (2020);pii:ciaa344. doi:10.1093/cid/ciaa344; 3. Wölfel R, et al. Nature (2020);[epub ahead of print]. doi:10.1038/s41586-020-291-x. 4. Jin Y, et al. Int J Infect Dis (2020). 5. Liu W, et al. J Clin Microbiol (2020);pii:JCM.00461-20.doi:10.1128/JCM.00461-20. 6. Guo L, et al. Clin Infect Dis (2020); doi:10.1093/cid/ciaa310. 7. Zhang W, et al. Emerg Microb Infect;9:386-389..

# SARS-CoV-2 infection

*Molecular and serology tests serve different use cases*



# Elecsys® Anti-SARS-CoV-2



## Intended use

Elecsys® Anti-SARS-CoV-2 is an immunoassay for the in vitro **qualitative detection of antibodies (including IgG)** to Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) in human serum and plasma. The test is intended as an aid in the determination of the immune reaction to SARS-CoV-2.

The electrochemiluminescence immunoassay “ECLIA” is intended for use on Elecsys® and **cobas e** immunoassay analyzers.

# Links to educational & technical material

*Our commitment to help put a stop to the COVID-19 pandemic*

## Real-time RT-PCR test

Detecting SARS-CoV-2 in currently infected patients with high sensitivity

The **cobas**<sup>®</sup> SARS-CoV-2 Test is a real-time RT-PCR test intended for the qualitative detection of SARS-CoV-2 in nasopharyngeal and oropharyngeal swab samples from patients.

This test can run on Roche's fully automated **cobas**<sup>®</sup> 6800 and **cobas**<sup>®</sup> 8800 Systems under Emergency Use Authorization. The test is also available for countries accepting the CE-mark.

[cobas<sup>®</sup> SARS-CoV-2 Test >](#)

[cobas<sup>®</sup> 6800 System >](#)

[cobas<sup>®</sup> 8800 System >](#)

## Antibody test

Detecting a patient's immune response to SARS CoV-2 with high specificity

The Elecsys<sup>®</sup> Anti-SARS-CoV-2 is an immunoassay for the in vitro qualitative detection of antibodies (including IgG) to SARS-CoV-2 in human serum and plasma. The test is intended as an aid in the determination of the immune reaction to SARS-CoV-2.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on **cobas e** immunoassay analyzers.

[Elecsys<sup>®</sup> Anti-SARS-CoV-2 >](#)

[cobas e 411 analyzer >](#)

[cobas e 601 module >](#)

[cobas e 602 module >](#)

[cobas e 801 module >](#)

<https://diagnostics.roche.com/global/en/c/covid-19-pandemic.html>

[https://dialog.roche.com/global/en\\_us/covid/run2.html](https://dialog.roche.com/global/en_us/covid/run2.html)


# Roche SARS-CoV-2 diagnostic testing



*ASLM echo session, July 10<sup>th</sup> 2020*

---


*Dr. Matthias Strobl, Senior Clinical Science Leader - Infectious Diseases*




# Serum antibody responses during infection by SARS-CoV-2

*Anti-SARS-CoV-2 IgM is not a specific indicator of early infection*


Antibodies are generated against the spike (S) and the nucleocapsid (N) protein <sup>1,2</sup>



IgM and IgG appear mostly at the same time, but IgG can rise earlier than IgM! <sup>2,3</sup>



Antibody levels are **not directly correlated** with clinical severity <sup>2-4</sup>



(1) Okba NMA, Müller MA, Li W, Wang C, et al. medRxiv 2020.03.18.20038059; (2) To KKW et al. (2020). Lancet Infect Dis 2020; DOI:10.1016/S1473-3099(20)30196-1 (3) Wu F et al. medRxiv 2020.03.30.20047365; doi: <https://doi.org/10.1101/2020.03.30.20047365>; (4) Long et al. (2020). Nat Med. <https://doi.org/10.1038/s41591-020-0897-1>

\*after symptom onset

# Nucleocapsid selected as target antigen

## *Elecsys<sup>®</sup> Anti-SARS-CoV-2*

### N selected based on experimental data and manufacturability:

- ✓ Performed better than all other targets in several hundred clinical samples
- ✓ Can be synthesized in soluble form
- ✓ Large-scale production in a fast and reliable manner

### Main performance requirement for selection:


- Y High specificity

**Nucleocapsid:**  
**Selected**

**Envelope:**  
4 prototypes

**Spike:**  
29 prototypes

**Membrane:**  
5 prototypes



# Considerations when testing for SARS-CoV-2 antibodies

## *Epidemiological screening*



For SARS-CoV-2 epidemiological screening, the aim is to **correctly identify** all patients who have been **exposed to the SARS-CoV-2** to better understand the prevalence, epidemiology and dynamics of the virus and COVID-19<sup>1,2</sup>

For this purpose, a **high positive predictive value** is highly desirable<sup>2</sup>



Given that **levels** and **chronological** order of IgM and IgG antibody appearance are **variable**, **targeting both antibodies simultaneously** following a total antibody approach may be beneficial for increasing the **sensitivity** of a test<sup>3,4</sup>

# Considerations when testing for SARS-CoV-2 antibodies<sup>1,2</sup>

## *Antibody maturity*

Ab Main Types →	Early/Immature	Mature	Neutralizing
<b>Description</b>	<ul style="list-style-type: none"> <li>Appear in <u>early infection</u> phase</li> <li><b>Do not effectively recognize the virus</b></li> </ul>	<ul style="list-style-type: none"> <li>Appear in <u>convalescent</u> phase</li> <li><b>Effectively recognize the virus</b></li> </ul>	<ul style="list-style-type: none"> <li>Appear in <u>convalescent/immunity</u> phase</li> <li><b>Effectively neutralize the virus</b></li> </ul>
<b>Examples</b>	IgM, early/immature IgA, IgG	Late/mature Igs	Neutralizing antibodies ( <i>sub-set of mature Igs</i> )
<b>Relevance/Purpose</b>	Initial host response to start understanding the virus	Host memory of the virus for future recognition	Render the virus ineffective against the host

**Note:** All neutralizing antibodies are mature antibodies BUT not all mature antibodies are neutralizing antibodies<sup>1,2</sup>


*Targeting mature antibodies may be beneficial in increasing testing specificity*

# Fit for Purpose Antibody Assay Formats




## Mature High Affinity Total Antibody Assay In-solution Double Antigen Sandwich<sup>2</sup>

Single arm (high affinity) binding  
is necessary to give positive result



Complementary antigens (Ag) are in solution

Early/immature Ab  
(No binding)




Assay requirement of 2 independent binding antibodies to complementary Ags enables specific, preferential detection of late/mature Abs

<sup>1</sup> Roche internal data; <sup>2</sup> Faatz, E., Finke, A et al. (2019). In *Analytical Electrogenerated Chemiluminescence* (pp. 443-470)

# SARS-CoV-2 neutralizing antibodies

## *Spike vs nucleocapsid as target antigen*

### How do antibodies inhibit viruses? <sup>1-4</sup>



### Neutralizing Abs may confer protection

Also considered critical for selection of convalescent plasma for therapeutic or prophylactic uses, and for vaccine assessment

Studies show **anti-nucleocapsid & anti-spike Ab** correlate with neutralization, i.e. putative clinical immunity <sup>9-11</sup>

**In the Elecsys<sup>®</sup> assay format, the nucleocapsid antigen confers high specificity**



*\* Examples of non-neutralizing complex interactions<sup>2-4</sup>: 1) Tagging cells for destruction by effector cells before viruses have a chance to exit (e.g. antibody-dependent cellular cytotoxicity (ADCC); 2) Induction of complement-mediated cytolysis, increased T cell responses associated with enhanced dendritic cell function, and reduced viral replication in culture; 3) Cooperatively elicit robust protective immunity when combined virus-specific CD8 T cells<sup>4</sup> NOTE: also cellular immunity plays a role<sup>6,7</sup>*

(1) Burton, D (2002). *Nat Rev Immunol* 2, 706–713; (2) Payne S (2017). *Viruses: Chapter 6 - Immunity and Resistance to Viruses*, Editor(s): Susan Payne, Academic Press; (3) Coughlin MM (2012). *Rev Med Virol*. 22(1):2–17.; (4) Zhou G et al. (2020). *Int J Biol Sci*;16(10):1718–1723; (5) Murin CD et al. (2019). *Nat Microbiol* 4, 734–747; (6) Damian M et al. (2008). *J Immunol*;181(6):4168–4176; (7) Laidlaw BJ et al. (2013). *PLoS Pathog* 9(3): e1003207; (8) Enayatkhani M et al. (2020). *PLoS Pathog* 9(3): e1003207; (9) To K et al. *Lancet Infect Dis*. 20(5), 565–74; (10) Wu F et al. (2020). *JAMA*. 323(13), 1239–1242; (11) Okba N et al. *medRxiv*. 2020. preprint doi: <https://doi.org/10.1101/2020.03.18.20038059>

# Correlation of Elecsys<sup>®</sup> Anti-SARS-CoV-2 with neutralization

*Positive correlation with neutralizing capacity*

- A positive test result with Elecsys Anti-SARS-CoV-2 correlates well with neutralization capacity in the same sample.
- Number of negative samples is too low to determine specificity, but no false positive results were observed with Elecsys<sup>®</sup> Anti-SARS-CoV-2.

**Roche internal study: 46 clinical samples from individual patients <sup>1</sup>**


		Pseudo-Neutralisation <sup>2</sup>		
		Pos	Neg	Total
Elecsys <sup>®</sup> Anti-SARS-CoV-2	Pos	38	0	38
	Neg	6	2	8
	Total	44	2	46
PPA		86.4 % (73.3 % - 93.6 %)*		
PNA		100.0 % (34.2 % - 100.0 %)*		
POA		87.0 % (74.3 % - 93.9 %)*		

PPA: percent positive agreement, PNA: percent negative agreement, POA: percent overall agreement; N: nucleocapsid, S: spike

\* 95% confidence interval; # antibody tests were performed 49 days after PCR.

1) Roche, data on internal file; 2) Meyer B, Torriani G, Yerly S, et al. Validation of a commercially available SARS-CoV-2 serological Immunoassay. medRxiv. 2020. <https://doi.org/10.1101/2020.05.02.20080879>

# What are the objectives for testing for antibodies to SARS-CoV-2



## What can antibody detection tell us?



Serologic assays that accurately assess prior infection to SARS-CoV-2 will be essential for **epidemiologic studies**, ongoing **surveillance**, and **vaccine studies**<sup>1</sup>



## What can't antibody detection tell us?




Negative results cannot exclude SARS-CoV-2 infection, particularly among those with recent exposure to the virus<sup>1</sup>



There is currently no conclusive evidence that people who have recovered from COVID-19 and have antibodies are protected from a second infection<sup>2</sup>


# When testing for antibodies.....

## *Importance of sensitivity*




### Sensitivity

**Highly sensitive** antibody tests:



High likelihood of correctly **identifying someone who was infected, as having been infected**



Avoid reporting someone who was infected, as having not been infected

# Elecsys® Anti-SARS-CoV-2

*High clinical sensitivity in patients with proven past infection (extended data set\*)*

**496** samples from:

---

**103** symptomatic patients PCR confirmed SARS-CoV-2 infection

---

Tested with **Elecsys® Anti-SARS-CoV-2**

One or more sequential specimens from these patients were collected after PCR confirmation at various time points

Days post PCR confirmation	N	Non-reactive	Sensitivity (95 % CI*)
0 – 6 days	161	64	60.2 % (52.3 – 67.8%)
7 – 13 days	150	22	85.3% (78.6 – 90.6%)
<b>≥ 14 days</b>	<b>185</b>	<b>1#</b>	<b>99.5% (97.1 – 100%)</b>

**Elecsys® Anti-SARS-CoV-2 has 99.5% sensitivity in detecting antibodies after the reported median day of seroconversion**





# Clinical sensitivity

## Public Health England assay evaluations

**Sensitivity of Elecsys® increases with time after symptom onset, indicating correlation with sero-conversion and/or increasing proportion of late, high-affinity antibodies.**

**Sensitivity of Architect at the same level across all time intervals indicating detection of low- and high-affinity antibodies.**

**Sensitivity of Vitros and Liaison is insufficient.**




\* unclear interval, as the date from admission into hospital was supplied rather than the date of symptom onset

1) Duggan J et al. (2020). 2) Duggane t al. (2020)


# When testing for antibodies.....

## *Importance of specificity*




 **Specificity**  
**Highly specific** antibody test:

---

 Avoid reporting someone who was infected, as having not been infected

---

 Avoid reporting someone who was not infected, as having been infected

# Elecsys® Anti-SARS-CoV-2



*Fit-for-purpose clinical specificity proven in a large sample cohort*

 **10533** samples from:

---

 Diagnostic routine

 Blood donors

 A common cold panel

 A coronavirus panel\* obtained before December 2019

---

 Tested with **Elecsys® Anti-SARS-CoV-2**

Cohort	N	Reactive	Specificity % (95 % CI)
Diagnostic routine	6305	12	99.81 % (99.67 – 99.90 %)
Blood donors	4148	9	99.78 % (99.59 – 99.90 %)
Common cold panel	40	0	100% (91.19 – 100%)
Coronavirus panel*	40	0	100% (91.19 – 100%)
<b>Overall</b>	<b>10453</b>	<b>21</b>	<b>99.80 % (99.69 – 99.88 %)</b>

**Elecsys® Anti-SARS-CoV-2 has 99.80% overall specificity.  
No cross-reactivity with the common cold coronaviruses was observed.**

\*40 potentially cross-reactive samples from individuals with past infection with coronavirus HKU1, NL63, 229E, or OC43, confirmed by PCR.  
Roche Diagnostics GmbH, R&D department; Penzberg, Germany. R&D record ACOV2-113-118\_sensi-speci Elecsys® Anti-SARS-CoV-2 method sheet (v1, Apr 2020).

# Elecsys® Anti-SARS-CoV-2



*Fit-for-purpose clinical specificity proven in a large sample cohort*



**747** potentially cross-reactive samples including:



**40 samples** from individuals with common cold symptoms



**40 samples** from individuals confirmed to be infected with one of the four common cold CoVs\*



Tested with **Elecsys® Anti-SARS-CoV-2**

Cohort	N	Reactive	Specificity % (95 % CI)
Common cold panel	40	0	100% (91.19 – 100%)
Coronavirus panel*	40	0	100% (91.19 – 100%)
Other potentially cross-reacting samples <sup>§</sup>	667	4 <sup>#</sup>	99.4% (95.6 – 99.8%)
<b>Overall</b>	<b>747</b>	<b>4</b>	<b>99.47% (98.63 – 99.85%)</b>


**No cross-reactivity with the endemic common cold coronaviruses was observed.**

\*from individuals with past infection with coronavirus HKU1, NL63, 229E, or OC43, confirmed by PCR; <sup>§</sup> Pre-pandemic samples with reactivity for various other indications, which could have an elevated potential for unspecific interference, were tested for reactivity in the Elecsys anti-SARS-CoV-2 assay; <sup>#</sup> acute CMV infection (IgM+, IgG+): 1; acute EBV infection (IgM+, VCA IgG+): 2; systemic lupus erythematosus: 1. Roche Diagnostics GmbH, R&D department; Penzberg, Germany. R&D record ACOV2-113-118\_sensi-speci

# Elecsys® Anti-SARS-CoV-2



*Only assay confirmed to not cross-react with endemic coronavirus samples*



	Roche	Abbott	DiaSorin	Biomérieux	OCD	Siemens	Beckman	Snibe	Mindray
<b>Test name</b>	Elecsys Anti-SARS-CoV-2	Architect SARS-CoV-2 IgG	Liaison SARS-CoV-2 S1/S2 IgG	Vidas SARS-COV-2 IgM	Vitros Anti-SARS-CoV-2 Total / IgG	Atellica SARS-CoV-2 Total	Access SARS-CoV-2 IgG	Maglumi 2019 n-CoV IgG	SARS-CoV-2 IgG
<b>HKU1, NL63, OC43 or 229E</b>	40 samples 	No data 	Only 4 samples 	No data 	No data 	No data 	No data 	No data 	Only 5 samples 
<b>Interference</b>	No interference 	Potential interference 	Missing 	Missing 	Potential interference 	Potential interference 	Potential interference 	Potential interference	Missing 

Elecsys Anti-SARS-CoV-2 Method sheet (April 2020); Snibe Maglumi 2019-nCoV IgG & IgM Method sheet (March 2020); Abbott Architect SARS-CoV-2 IgG Method Sheet (April 2020)  
DiaSorin LIAISON® SARS-CoV-2 S1/S2 IgG Method Sheet (April 2020); Mindray SARS-CoV-2 IgG & IgM (CLIA) Method Sheet (March 2020); Ortho Clinical Diagnostics VITROS Immunodiagnostic Products Anti-SARS-CoV-2 Total Reagent Pack Method Sheet (April 2020); Shenzhen Yhlo Biotech iFlash-SARS-CoV-2 IgG & IgM Flyer; Qian Ch, Zhou M, et al. 2020. Development and Multicenter Performance Evaluation of The First Fully Automated SARS-CoV-2 IgM and IgG Immunoassays. medRxiv preprint doi: <https://doi.org/10.1101/2020.04.16.20067231>.

# Head-to-head comparison of three automated antibody assays

*Elecsys<sup>®</sup> with highest PPV at low prevalence*

- While all three assays presented with high specificities, at low prevalence, the minor differences in specificity resulted in profound discrepancies of positive predictability.
- At 1% prevalence, the PPV values of Roche and DiaSorin differ so clearly that not even the 95% CI intervals overlap.





# The importance of a high-specificity test

*A positive test result in New York is different from a positive test result in Santa Clara...*



A test's **specificity** plays an increasingly important role in populations where the infection rate (or prevalence) **is lower**.<sup>1,2</sup>



\*The positive predictive value is the probability that persons with a positive test result truly have the disease.

1. Bonislowski, A. (2020, April 27). New York, California Serology Studies Give Early Estimates of COVID-19 Prevalence. Retrieved from: [https://www.360dx.com/infectious-disease/new-york-california-serology-studies-give-early-estimates-covid-19-prevalence?utm\\_source=Sailthru&utm\\_medium=email&utm\\_campaign=360DN%20Tues%202020-04-28&utm\\_term=360Dx%20Daily%20News#.Xqv\\_4qgzZaR](https://www.360dx.com/infectious-disease/new-york-california-serology-studies-give-early-estimates-covid-19-prevalence?utm_source=Sailthru&utm_medium=email&utm_campaign=360DN%20Tues%202020-04-28&utm_term=360Dx%20Daily%20News#.Xqv_4qgzZaR)

2. Saplakoglu, Y., Writer, S. (2020, April 23). 1 in 5 people tested in New York City had antibodies for the coronavirus. Retrieved from: <https://www.livescience.com/covid-antibody-test-results-new-york-test.html>

3. Bendavid, E. et al. (2020, April 11). COVID-19 Antibody Seroprevalence in Santa Clara County, California. medRxiv preprint doi: <https://doi.org/10.1101/2020.04.14.20062463>

# The importance of a high-specificity test

*A positive test result in New York is different from a positive test result in Santa Clara...*



A test's **specificity** plays an increasingly important role in populations where the infection rate (or prevalence) is **lower**.<sup>1,2</sup>



## New York

100,000  
people tested

15,000  
infected with  
SARS-CoV-2

15%  
prevalence<sup>1,2</sup>



Other on-market antibody  
test with **98.80%**  
specificity




Elecsys<sup>®</sup> Anti-SARS-CoV-2  
with **99.80% specificity**



16,020 positive results



15,170 positive results



**6%** probability that the positive  
result is **incorrect**  
**94%** probability that the  
positive result is **correct\***



**1%** probability that the positive  
result is **incorrect**  
**99%** probability that the positive  
result is **correct\***



## Santa Clara

100,000  
people tested

1,500  
infected with  
SARS-CoV-2

1.5%  
prevalence<sup>3</sup>



Other on-market  
antibody test with **98.80%**  
specificity



Elecsys<sup>®</sup> Anti-SARS-CoV-2  
with **99.80% specificity**



2,682 positive results



1,697 positive results



**44%** probability that the  
positive result is **incorrect**  
**56%** probability that the  
positive result is **correct\***



**12%** probability that the positive  
result is **incorrect**  
**88%** probability that the positive  
result is **correct\***

\*The positive predictive value is the probability that persons with a positive test result truly have the disease.

1. Bonislowski, A. (2020, April 27). New York, California Serology Studies Give Early Estimates of COVID-19 Prevalence. Retrieved from: [https://www.360dx.com/infectious-disease/new-york-california-serology-studies-give-early-estimates-covid-19-prevalence?utm\\_source=Salithru&utm\\_medium=email&utm\\_campaign=360DN%20Tues%202020-04-28&utm\\_term=360Dx%20Daily%20News#.Xqv\\_4qgzZaR](https://www.360dx.com/infectious-disease/new-york-california-serology-studies-give-early-estimates-covid-19-prevalence?utm_source=Salithru&utm_medium=email&utm_campaign=360DN%20Tues%202020-04-28&utm_term=360Dx%20Daily%20News#.Xqv_4qgzZaR)

2. Saplakoglu, Y., Writer, S. (2020, April 23). 1 in 5 people tested in New York City had antibodies for the coronavirus. Retrieved from: <https://www.livescience.com/covid-antibody-test-results-new-york-test.html>


3. Bendavid, E. et al. (2020, April 11). COVID-19 Antibody Seroprevalence in Santa Clara County, California. medRxiv preprint doi: <https://doi.org/10.1101/2020.04.14.20062463>

# How does Elecsys® Anti-SARS-CoV-2 serology test provide value

## *An antibody detection test*



The **Elecsys®** Anti-SARS-CoV-2 **immunoassay** is an *in vitro* test, using human serum and plasma drawn from a blood sample, to **detect total antibodies** (IgA, IgM and IgG) indicating an immune response to SARS-CoV-2



The test may be used:



In **epidemiological research** to help better understand the spread of the disease, especially in people who may have been **infected but did not display symptoms**



Together with **molecular tests** to aid in the *post hoc* diagnosis of suspected COVID-19 patients



Hospitals and reference laboratories can run the test on Roche's **cobas e** analysers, which are **widely available in laboratories around the world**

# What are the unanswered questions for COVID-19 / SARS-CoV-2?

*Do antibodies confer reliable immunity?*

*How long do antibodies confer immunity?*

*How long is a patient infective?*

*From when on is a patient protected?*

*Doing now what patients need  
next*