

G
A
D Humanity
Tested

What does COVIOS Ag offer?

Introducing a new, high quality COVID-19 antigen test from Global Access Diagnostics

Mark Radford, Executive Director at Global Access Diagnostics (GAD)

Dr. Emily Adams, Director of Epidemic and Neglected Diseases at Mologic & GAD

September 2021

Contents

- ✓ **Company introduction**
- ✓ **COVIOS Rapid Antigen Test product overview**
- ✓ **Clinical studies & performance**
- ✓ **Regulatory status**
- ✓ **Product pipeline**

Company Introduction

G
A
D Humanity
Tested

COVID-19 has exposed the difficulties lower- and middle-income countries (LMICs) face in obtaining rapid diagnostic tests (RDTs). With global demand outstripping supply, poorer countries often find themselves at the back of the queue. GAD was founded to address this issue.

We are...

- a **social enterprise**, spun out of Mologic in 2020 with the goal of **building access to reliable, affordable RDTs in LMICs**.
- funded through the **Bill and Melinda Gates Foundation, Foundation for Innovative New Diagnostics** and **Soros Economic Development Fund**. Longer term, GAD will be wholly owned by the charity 'Global Access Health'.
- **Re-invested profits** generated back into the service of our mission
- granted **perpetual, royalty-free access to IP** generated by Mologic to develop improved RDT platforms and tests for neglected endemic and epidemic infections for LMICs.
- expanding **manufacturing capacity** to make **>200 million RDTs** by the end of **2021**.
- working to **support the development of diagnostic manufacturing capacity directly within LMICs** and build financially sustainable entities that are not profit driven and responsive to local needs.
- a **developer-agnostic manufacturing platform**. We aim to become an active partner to the donor, philanthropic and LMIC community to help address RDT innovation and supply gaps.

See *Financial Times* coverage [here](#)

To find out more, contact us on info@globalaccessdiagnostics.com

G
A
D **Humanity**
Tested

The world's largest social enterprise for diagnostics



**OPEN SOCIETY
FOUNDATIONS**

Confidential

History of Global Access Diagnostics

- 1988** Clearblue pregnancy test developed as first commercial lateral flow device (LFD). Prof. Paul Davis named one of three inventors.
- 2003** Paul Davis his son Mark established Mologic, an R&D company focused on LFD technology.
- 2016** Investment from BMGF establishes CARD laboratory at Mologic, dedicated to innovation for epidemic and neglected diseases.
- 2020** Global Access Diagnostics (GAD) spun out of Mologic as a social enterprise LFD manufacturing platform focused on the needs of low and middle-income countries. GAD is funded by BMGF, FIND, Soros Economic Development Fund (SEDF) at the Open Society Foundations, and UK government.
- 2021** Mologic was bought by Global Access Health (GAH). The deal was funded by a group of impact investors led by SEDF. GAH is a not-for-profit holding company that also owns GAD.

Our North Star:

The best tests to the most people at the lowest cost



**OPEN SOCIETY
FOUNDATIONS**

What differentiates Global Access Diagnostics?

1. We are developer agnostic & LMIC focused

We have a strong internal product pipeline, but we are also happy to manufacture tests for other developers, when these meet our objectives of quality, relevance and affordability for low and middle-income countries (LMICs).

2. We are committed to transparent pricing

For LMICs, we operate on a COGS+ basis, meaning that our objective is to only cover costs and a small margin for sustainability, rather than maximizing profit.

3. We support distributed manufacture

We are working with partners in several LMICs to move the center of gravity for diagnostic test manufacturing to those regions.

4. We believe in innovation

We have projects ongoing not only to innovate in the science, but also in areas such as LFD manufacturing technology, eco-appropriate materials and data enablement.

COVIOS Rapid Antigen Test product overview

$\frac{\text{G}}{\text{A}} \frac{\text{D}}{\text{D}}$ Humanity
Tested

Product Features of GAD COVIOS Ag

COVID 19 RAPID ANTIGEN TEST

**Global
Access
Diagnostics**

Feature	Details
Brand name	<ul style="list-style-type: none"> COVIOS
Intended use	<ul style="list-style-type: none"> Point of care LFA for qualitative detection of SARS-CoV-2 Ag Individuals with symptoms during: <ul style="list-style-type: none"> Acute phase of infection (no cap on days included) Pre-symptomatic or asymptomatic but suspected of COVID-19 by provider
Clinical data	<ul style="list-style-type: none"> Sensitivity: 90.6% (85.6% to 94.0%) Specificity: 100.0% (99.2% to 100.0%)
Detection of variants	<ul style="list-style-type: none"> Reliable detection of the following variants: <ul style="list-style-type: none"> B1.1.7 in UK, B1.351 in South Africa, and P1 in Brazil
Limit of detection	<ul style="list-style-type: none"> 2.5×10^2 pfu/ml
Sample type	<ul style="list-style-type: none"> Nasal swabs
Conditions	<ul style="list-style-type: none"> Stable storage at room temperature (2-30°C)
Shelf-life	<ul style="list-style-type: none"> 18 months (24 months expected by September 2021)
Time to results	<ul style="list-style-type: none"> 10 minutes
Pack size	<ul style="list-style-type: none"> 25 tests per pack (individually wrapped)
Pack contents	<ul style="list-style-type: none"> Test kit: <ul style="list-style-type: none"> Lateral flow device Swab extraction buffer tube and nozzle Sterile swab Instructions for use Positive control



Confidential

User guide: Specimen Collection and Handling

1 Open kit components pictured above and lay on a clean, flat surface

Sterile swab

Lateral flow test device

Sample hole T C Read window

Swab Extraction Buffer Tube

2 NASAL SWAB SPECIMEN COLLECTION:

1. Ask patient to blow their nose thoroughly to remove excess mucus.
2. Insert swab 1-2cm into both nostrils one at a time and rotate the swab against the nasal wall 5-6 times and withdraw.

X5-6

3 TUBE: Peel back and remove the foil, Lower the swab specimen into the tube. Submerge the swab in the buffer and mix thoroughly pressing the swab on the tube sides for 5 to 10 seconds.

5-10 SECONDS

4 SQUEEZE: Extract as much fluid as possible by squeezing and twisting the swab head through your fingers as you draw it out from the tube.

PRESS ON THE SIDE OF THE TUBE

CLOSE: Seal the tube by pushing the nozzle in place.

5 Ensure tube is **VERTICAL** and carefully dispense 5 continuous drops of extract into the sample hole on the lateral flow test device by gently squeezing the tube.

x5

10 min

⚠ Read the TEST (T) visually at 10 minutes.

⚠ Discard extraction tube and swab.

INTERPRETATION OF TEST RESULTS

⚠ DO NOT INTERPRET THE TEST AFTER 15 MINUTES. IF UNSURE, PLEASE REPEAT TEST ON NEW DEVICE AND INTERPRET AT 10 MINUTES

POSTIVE	NEGATIVE	INVALID
<p>T C</p> <p>T C</p> <p>If the test line (T) is visible, this indicates a positive SARS-CoV-2 result.</p> <p>If the test line (T) is faint, this indicates a positive SARS-CoV-2 result.</p>	<p>T C</p> <p>If the test line (T) is not visible, this indicates a negative result.</p>	<p>T C</p> <p>If the control line (C) is absent, the test is invalid.</p>

Clinical studies & performance

$\frac{\text{G}}{\text{A}}$
 $\frac{\text{A}}{\text{D}}$

Humanity
Tested

Global clinical sites for the validation of Mologic/GAD COVID-19 Rapid Antigen Test POC Cassette



FIND/WHO

Dr Claudia Denkinger (Germany)
UHCP (Peru)

Institut Pasteur de Dakar, Senegal

Dr Amadou Sall
Dr Cheikh Tidiane Diagne

Zankli Research Centre, Nigeria

Dr John Bimba
Prof Luis Cuevas

University of Witwaterstand, South Africa

Prof Wendy Stevens
Prof Lesley Scott

Instituto Nacional de Saude (INS), Mozambique

Dr Ilesh Jani

Federal University of Sergipe, Brazil

Prof. Ricardo Gurgel

UK

Dr Emily Adams, LSTM Dr
Nicholas Easom, HUTH Dr
David Tate, NHNFT
Dr Rahul Batra, GSTT
Dr Tim Planche, SGUL
Prof. Sanjeev Krishna, SGUL

USA

Dr Paul Drain, UW
Dr Tyler Miller, MGH

Europe:

UK
Germany

East Asia:

Indonesia
Malaysia

Americas:

USA
Peru

Africa:

Senegal
Nigeria
South Africa
Mozambique



Funders



GAD COVIOS Ag Rapid Test: Clinical results of FIND independent evaluation

- GAD COVIOS Ag rapid diagnostic test (developed at Mologic) has demonstrated comparable performance to lateral flow with WHO emergency use listing in independent evaluations by FIND
- GAD cassette format have been submitted to the WHO for EUL and USFDA for EUA

Lateral flow tests with WHO EUL

	Location	Setting	% of patients symptomatic	Total N	Positivity	Specificity	Sensitivity by cycle threshold value					
							N	All	N	≤33	N	≤25
	Germany	Drive-in	46.2%	281	16% (44/281)	99.2%	44	86.4% (73.3-93.9%)	42	90.5% (77.9-96.2%)	31	96.8% (83.8%-99.4%)
	Germany	Ambulatory	96.6%	179	23% (41/179)	99.3%	41	80.5% (66-89.8%)	32	87.5% (71.9-95%)	21	100% (84.5-100%)
	Brazil	Community	100%	214	36% (78/214)	99.3%	78	84.6% (75-91%)	72	91.7% (83-96.1%)	46	100% (92.3-100%)
	Germany	Drive-in + Ambulatory	63.1%	529	19% (100/529)	97%	100	91% (83.8-95.2%)	96	93.8% (87-97.1%)	80	97.5% (91.3-99.3%)
	India	Hospital	20.8%	600	18% (105-600)	99.6%	105	74.3% (65.2-81.7%)	103	74.8% (65.6 - 82.2%)	81	87.7% (78.7-93.2%)
	Germany	Ambulatory and drive-in	66.5%	665	29% (194/665)	100%	191	90.6% (85.6-94%)	186	92.5% (87.8-95.5%)	166	96.4% (92.3-98.3%)

Confidential

Source: FIND Evaluation External Reports

Top performance in FIND Germany trial for WHO EUL

FIND independent evaluation of 665 participants in Germany found a 91% sensitivity and 100% specificity

Information from IFU (Source: Germany FIND study)	Prospective Recruitment		RT-qPCR		
			Pos	Neg	Total
	GAD Rapid Antigen Test	Pos	173	0	173
		Neg	18	458	476
Total		191	458	649	
Performance	Sensitivity		Specificity		
	90.6% (85.6% to 94.0%)		100.0% (99.2% to 100.0%)		
<ul style="list-style-type: none"> • GAD COVID-19 Rapid Antigen Test uses a nose only (anterior nares) swab • RT-qPCR platforms: TibMolbiol and Roche • Range Ct: 11.7 to 34.7 • FINDDX Ref: https://bit.ly/3nEei2 					

Sensitivity stratified by RT-PCR cycle threshold value	
Ct	Sensitivity
<20	100.0% (92/92)
<25	96.4% (160/166)
<33	92.5% (172/186)

Information directly from FIND report	Lowest dilution detected	Verified LOD concentration	Viral Copy equivalent	Supplier-reported LOD	Reference PCR method	<ul style="list-style-type: none"> • LightMix® Modular SARS-CoV (COVID19) E-gene (Tib Molbiol) <ul style="list-style-type: none"> ○ N = 323 • Cobas SARS-CoV-2 (Roche Diagnostics Inc) <ul style="list-style-type: none"> ○ N = 342
	Analytical Sensitivity	2.5 x10² pfu/ml ~ 3.52 x 10² TCID₅₀/ml	2.5 x10² pfu/ml	5.9 x10⁵ copies/ml applied to test		

Note: viral dilution was applied directly to the test cassette, not to the provided swab

Regulatory status

$\frac{\mathbf{G}}{\mathbf{A}} \frac{\mathbf{D}}{\mathbf{D}}$ Humanity
Tested

Regulatory Strategy for GAD COVID-19 Rapid Antigen Test

Aggregate regulatory pathways

- **CE Mark**
 - ✓ Achieved in December 2020
- **MHRA (UK)**
 - ✓ Achieved in December 2020
- **WHO EUL**
 - ✓ Submitted in May 2021: [link](#)
 - Responding to first round of WHO questions
- **FDA EUA**
 - ✓ Received RADx funding in April 2021
 - ✓ Submitted application in June 2021
 - Awaiting approval


In-country validation & verification

Region	Country	V&V organizations
Africa	Botswana	Ministry of Tertiary Education, Research, Science and Technology (contact from President's Office)
	Burkina Faso	Institut de Recherche En Sciences de la Sante (IRSS)
	Cameroon	National Research Lab (NRL) or Centre International de Référence Chantal Biya (CIRCB) or Centre Pasteur du Cameroun (CPC)
	DRC	Institut National pour la Recherche Biomedicale (INRB)
	Ethiopia	Ethiopian Public Health Institute (EPHI)
	Ghana	Noguchi Labs/ Public Reference lab
	Kenya	Kenya Medical Research Institute (KEMRI)
	Lesotho	National Research Lab (NRL)
	Malawi	Department of HIV and AIDS, Ministry of Health (CHSU) - (National Research Lab)
	Mali	Institut National de Santé Publique (INSP)
	Mozambique	Instituto Nacional de Saúde (INS) / Liverpool School of Medicine
	Nigeria	Medical Laboratory Science Council of Nigeria (MLSCN); Zankli Research Centre with funding from Wellcome Trust (<i>V&V not needed</i>)
	Rwanda	Rwanda Biomedical Centre (National Research Lab - NRL)
	Senegal	Institut Pasteur de Dakar
	South Africa	University of Witwatersrand
	LATAM	Brazil
Colombia		Instituto Nacional de Salud
Peru		FIND
Rest of World	Bangladesh	Barisal Biotech
	India	Indian Council of Medical Research (ICMR)
	Indonesia	University of Indonesia
	Pakistan	Aga Khan University
	Philippines	Biological Life Science
	Timor-Leste	<i>V&V not needed</i>

Confidential

Product pipeline

$\frac{\mathbf{G}}{\mathbf{A}} \frac{\mathbf{D}}{\mathbf{D}}$ Humanity
Tested



Contact details

info@globalaccessdiagnostics.com

G
—
A
—
D Humanity
Tested