



Innovative solutions for TB infection screening in decentralised settings

Davide Manissero, MD, MRCPCH, MSc, DTM&H

QIAGEN, Chief Medical Officer


VP, Head of Clinical, Medical, Regulatory Affairs , Medical Affairs

Diagnostic are essential for universal health coverage




47%

Of people around the world still don't have access to basic diagnostics



Current TB infection care cascade gaps



- **Steps in the cascade associated with greater losses included completion of testing**
- Among the reasons associated with losses for completing screening and testing:
 - Health systems issues (lack of human resources, hard to access clinic, health insurance etc..)
 - Social issues (wrong perception of the disease, stigma, mistrust etc...)

Implementation of TB infection Test & Treat is essential for TB elimination

TODAY



Care site



Patients empirically treated



Acceptability
Test accessibility



TEST



TREAT

Care site



TOMORROW



TEST



TREAT

Care site



Acceptability
Test accessibility

We want to bring the potential of QFT-plus performances in low- and middle-income countries and make it **accessible** to those in need

QIAreach QuantiFERON-TB



Evolution of QFT Technology

First generation
QuantiFERON-TB



2001: FDA approved

- Measured cell-mediated immunity to tuberculin purified protein derivative (PPD)
- **Breakthrough: TST becomes a blood test**

Second generation
QuantiFERON-TB
Gold



2004: FDA approved

- Liquid antigen version
- Antigens specific for M.tb with 99% specificity
- **Clinical benchmark: No cross reactivity with BCG**

Third generation
QuantiFERON-TB
Gold In-Tube



2007: FDA approved

- Logistical advantage – remote incubation
- **Lab benchmark: Scalable**
- >1500 peer reviewed publications
- >30 million tests sold

Fourth generation
QuantiFERON-TB
Gold Plus (QFT-Plus)



2014: CE-IVD
2017: FDA approved

- **Addition of patented CD8 antigens** – potential biomarker of intracellular TB burden
- New flexible blood draw options

Latest generation offer

QIAreach
QuantiFERON-TB*



Expected in Q3 2021

- **Using QFT-Plus technology to increase access to IGRA testing in high burden / Low resource settings**

* Currently under development

QIAreach Solution to address testing needs in low resource settings

No maintenance required


- No maintenance or calibration of the eHub
- Dust protection guard provided and easy to clean ports

Scalable solution

- Scalable from 1 to 8 samples, up to 24 samples/h

Up to 8 h battery backup

Allow uninterrupted operation when the power goes out or an external power source is not available



Digital and connected

- Optional software for results traceability
- Ability to send the results to LIMS

No cold chain

- No cold chain required for any test components

Portable

- eHub < 1 kg and portable

BLOOD COLLECTION

SEND TO LAB / STIMULATION

IGRA TEST


RESULTS & INTERPRETATION

Convenient options

4th generation T cell stimulation technology

TB infection testing made easy

High quality results





Clinical performances (Japan)

Objective

- To compare the QIAreach QuantiFERON-TB (QIAreach QFT) vs. QuantiFERON-TB Gold Plus assay (QFT-Plus) to detect tuberculosis (TB) infection;
- To evaluate diagnostic sensitivity of QIAreach QFT using active TB as surrogate for TB infection;
- To preliminarily evaluate QIAreach QFT in immunocompromised individuals

Results

	Healthy controls		Active TB		Total
	QFT-Plus positive	QFT-Plus negative	QFT-Plus positive	QFT-Plus negative	
Positive QIAreach QFT result	0	1	41	0	42
Negative QIAreach QFT result	0	41	0	0	41
Total	0	42	41	0	83



Conclusion

- This study demonstrates that QIAreach QFT test has high clinical performance: 100% sensitivity, 97.6% specificity, and **98.8% overall concordance using QFT-Plus as the reference standard**
- There is a statistically significant relationship between levels of IFN- γ in plasma of active TB patients and TTR suggesting that TTR could be used as a surrogate marker of IFN- γ concentration in plasma when using QIAreach QFT assay
- Seven cases in the active TB group who were immunocompromised (CD4 <200/ μ L) returned positive results on QIAreach QFT

Clinical performances (Internal study)

Objective

Agreement between QuantiFERON-TB Gold Plus and QIAreach QuantiFERON-TB: Performance of QIAreach QuantiFERON-TB (QIAreach QFT) was compared to QuantiFERON-TB Gold Plus (QFT-Plus) in a population with a mix of risk factors for TB infection. Specimens were collected from a total of 4 sites. All QFT-Plus ELISA testing and QIAreach QFT testing was performed at a single site. A total of 225 samples were included in the final performance comparison


Preliminary performance assessment

Counts	Frequency	Agreement	Upper 95% CI	Lower 95% CI
OPA*	219/225	97.3%	99.0%	94.3%
PPA	71/75	94.7%	98.5%	86.9%
NPA	148/150	98.7%	99.8%	95.3%

OPA: Overall percent agreement; PPA: Positive percent agreement; NPA: Negative percent agreement

* When factoring in 15 QFT-Plus indeterminate results, the OPA between QFT-Plus and QIAreach QFT is 91.3% (95% CI: 86.9 – 94.5%).

Plot of TB1-Nil and TB2-Nil values against QIAreach QFT Time for Positive Results



Conclusion

- In a preliminary performance evaluation versus QFT-Plus (n = 225), QIAreach QFT showed an OPA of 97.3% (95% CI: 99.0 – 94.3%), with a PPA of 94.7% (95% CI: 98.5–86.9%) and an NPA of 98.7% (95% CI: 99.8% – 95.3%).
- Significant correlation was observed between the QIAreach QFT time to result and IU/ml responses of QFT-Plus-positive samples.
- **QIAreach QFT shows a high level of agreement with QFT-Plus** and has the potential to overcome key hurdles for TB screening in high-burden, low-resource settings.


Usability study (Zambia)

Objective

Assess usability the QIAreach QFT and training needs for the assay implementation were assessed across three domains: (1) effectiveness, (2) efficiency, and (3) user satisfaction



Preliminary performance assessment



Effectiveness
Pass rate of 4 tasks



Efficiency
Completion of all tasks in less than 1 hour



Satisfaction
Likert scale 1-5



Conclusion

The characteristics of the platform together with our usability finding make the QIAreach-QFT assay suitable to be implemented in the remote area where limited infrastructure has hampered the accessibility of IGRA technologies to those in needs.

QFT-Plus demonstrated comparative cost-effectiveness with skin-based tests; high potential for QIAreach to exceed

Article | [Open Access](#) | Published: 11 December 2020

Cost-effectiveness of newer technologies for the diagnosis of *Mycobacterium tuberculosis* infection in Brazilian people living with HIV

Ricardo E. Steffen [✉](#), Marcia Pinto, Afranio Kritski & Anete Trajman

Scientific Reports **10**, Article number: 21823 (2020) | [Cite this article](#)

Strategy	Cost	Incremental cost	QALY	Incremental QALY
Diaskintest	884.70		8.386	
EC skin test	886.60	1.90	8.386	0
QFT-Plus	902.10	17.40	8.385	- 0.00055
TST PPD RT 23	925.50	40.80	8.356	- 0.02967

- Some operational aspects can have great impact on final test costs and, consequently, its cost-effectiveness
- No societal costs were included (indirect costs, loss of productivity or cost of death)
- No cost-effectiveness of opportunity cost were included

QIAreach – QFT: Hierarchy of effectiveness



