


SARS-CoV-2 Variants Tracking using Thermo Fisher Scientific qPCR Technology

Tshego Mothatinyane

Technical Sales Specialist

16/09/2021

 The world leader in serving science



Agenda

1 Background: Potential Variant Impact

2 Introducing: SNPs detection using TaqMan SARS-CoV-2 Mutation Panel

3 Features and qPCR Workflow


4 Technical Data

5 Questions


SARS-CoV-2 Viral Mutations

Viruses mutate. RNA viruses, like SARS-CoV-2, mutate at high rates in response to selective pressures


Continued uncontrolled transmission of SARS-CoV-2 in many parts of the world is creating conditions for significant virus evolution



SARS-CoV-2 has been mutating at a rate of about one to two mutations per month*



Some recently identified variants, however, have acquired mutations much more rapidly than scientists expected



*Nature <https://www.nature.com/articles/d41586-020-02544-6>

Potential Implications of New SARS-CoV-2 Variants

Potential Variant Impacts



Speed of
human-to-
human
transmission



Change
disease
severity



Susceptibility
to therapeutic
agents
(i.e., monoclonal
antibodies)



Evade
vaccine-
induced
immunity



Impact on
detection tests

SNPs detection in the current multi-variant analysis

Overview

- SNPs – Single Nucleotide polymorphisms detection of gene variations in a well.
- RUO Reflex technology; samples initially processed using TaqPath COVID-19 CE-IVD assay kits can be reflex tested using the same workflow. (LoD: Ct value of 30 or less)
- If the sample is positive and the S-gene is undetectable (S-gene “dropout”) the sample can be reflex-tested to RT PCR
- Multiple primer/ probe sets available to confirm mutations associated with 69/70del and B.1.1.7 or B.1.351, B.1.617.2 strains
- GenoTyping software using qPCR platforms (Any real-time PCR instrument such as the Applied Biosystems 7500, 7500 Fast, 7500 Fast Dx, QuantStudio5 and QuantStudio7 systems)

Global Surveillance of SARS-CoV-2 Variants

Multiple surveillance solutions to identify, detect, and confirm new and emerging SARS-CoV-2 variants and strains.

Public health partners



- Monitor viral genome
- Detect & Inform on emerging mutations
- Determine impact on detection, spread and vaccine / therapy effectivity


Genetic Surveillance

Next Gen Sequencing (NGS)

- Sequence full viral genome
- Detect new and emerging mutations

Capillary Electrophoresis

- Sequence targeted genes / areas of interest on viral genome (i.e., S-gene)



Mutation Verification

Capillary Electrophoresis

- Sequence verification of unexpected results
- Confirmation of mutations, strains or lineages

RT-qPCR

- Convenience of using same instrument for detection and mutation confirmation
- Customizable solution with a menu of assays

RT qPCR Detection

RT-qPCR

- Assess emerging mutations impact on performance claims
- Informs product development
- Drives rapid, proactive product modifications and communications


Research

Clinical

TaqMan® SARS-CoV-2 Mutation Panel

Convenient, customizable **real-time PCR** research solution for detection analysis and surveillance

- **Customizable panel**—build your own custom panel from a menu of verified assays to identify currently relevant SARS-CoV-2 mutations and adapt quickly as others emerge
- **Convenient**—use your current real-time PCR instrument to conduct follow up testing of SARS-CoV-2 samples
- **Scalable**—run a few or hundreds of samples to identify for one or many mutations
- **Unique, streamlined workflow**—combining our gold-standard TaqMan SNP Genotyping Assays with a 1-step RT-PCR reaction, go from RNA to results in just over 1 hour



TaqMan® SARS-CoV-2 Mutation Panel

PRODUCT DETAILS

Targets	<ul style="list-style-type: none">• Choose from a menu of verified real-time PCR SNP assays to build your custom panel
Assay design	<ul style="list-style-type: none">• Each assay includes two TaqMan minor groove binder (MGB) probes with nonfluorescent quenchers (NFQ):<ul style="list-style-type: none">• One VIC™ dye-labeled probe to detect the reference sequence• One FAM™ dye-labeled probe to detect the mutation sequence
Available sizes	<ul style="list-style-type: none">• 374 reactions• 1,250 reactions
Format	<ul style="list-style-type: none">• 1 mutation assay per well
Sample input	<ul style="list-style-type: none">• RNA extracted from SARS-CoV-2 samples with a CT value of less or equal to 30
Turnaround time	<ul style="list-style-type: none">• 1 hour and 10 min from extracted RNA to results
Recommended instruments	<ul style="list-style-type: none">• Any real-time PCR instrument such as the Applied Biosystems 7500, 7500 Fast, 7500 Fast Dx, QuantStudio5 and QuantStudio7 systems.
Recommended analysis software	<ul style="list-style-type: none">• QuantStudio™ Design and Analysis Software v2.5 or later with the Genotyping Analysis Module
Optional controls	<ul style="list-style-type: none">• AcroMetrix™ Coronavirus 2019 RNA Control• GeneArt™ DNA string or plasmid controls

Build a custom panel by choosing from our menu of over 50 verified SNP assays

Mutation	Gene	Associated Variants	MTO/MTS
A1708D	ORF1	B.1.1.7	MTO
A222V	S	B.1.177	MTO
A570D	S	B.1.1.7	MTO
A701V	S	B.1.351	MTO
D215G	S	B.1.351	MTO MTS
D614G	S	B.1.1.207, P.1, B.1.1.33, B.1.1.7, B.1.177, B.1.258, B.1.351, B.1.525, Mink Variant	MTS
D80A	S	B.1.351	MTO
delH69V70	S	B.1.1.7, B.1.258, B.1.525	MTS
delY144	S	B.1.1.7	MTS
E484K	S	P.1, B.1.1.33, B.1.351, B.1.525	MTS
E484Q	S	B.1.617.1, B.1.617.3	MTS
K417N	S	B.1.351	MTS
K417T	S	P.1	MTS

Mutation	Gene	Associated Variants	MTO/MTS
L18F	S	P.1, B.1.351	MTO MTS
L242_244L	S	B.1.351	MTO
L452R	S	B.1.617, B.1.429	MTS
N439K	S	B.1.258	MTO MTS
N501Y	S	P.1, B.1.1.7, B.1.351	MTS
P681H	S	B.1.1.207, B.1.1.7	MTS
P681R	S	B.1.617	MTS
Q27stop	ORF8	B.1.1.7	MTO
R246I	S	B.1.351	MTO
S982A	S	B.1.1.7	MTO
T20N	S	P.1	MTS
T716I	S	B.1.1.7	MTO
Y453F	S	Mink Variant	MTO

Key mutations and their distribution

SARS-CoV-2 | Variants and mutations

ThermoFisher
SCIENTIFIC

As SARS-CoV-2 continues to mutate, there are numerous variants popping up around the world. As new variants arise, stay up to date through our **Notable Variants** blog.

Did you know you can identify each of these key mutations using our TaqMan SARS-CoV-2 Mutation Panel? [Create your custom panel](#)



Alpha (B.1.1.7) B.1.177, B.1.258

First detected: United Kingdom

Key mutations:

- D614G Thought to make the coronavirus more infectious
- delH69V70 Alters the shape of the spike and may help it evade some antibodies
- delY144 Alters the shape of the spike and may help it evade some antibodies
- N439K Impacts ability to evade antibody-mediated immunity
- N501Y Helps the virus latch on more tightly to human cells
- P681H Predicted to enhance systemic infection and associated with increased transmissibility



Beta (B.1.351) B.1.1.33

First detected: South Africa

Key mutations:

- D215G May help the virus latch on tighter
- D614G Thought to make the coronavirus more infectious
- E484K Reduces antibody recognition; Associated with vaccine resistance
- K417N Helps the virus bind more tightly to human cells
- L18F May help the virus latch on tighter
- N501Y Helps the virus latch on more tightly to human cells



Gamma (P.1)

First detected: Brazil (and in Japan, detected in travelers from Brazil)

Key mutations:

- D614G Thought to make the coronavirus more infectious
- E484K Reduces antibody recognition; Associated with vaccine resistance
- K417T May help the virus latch on tighter
- L18F May help the virus latch on tighter
- N501Y Helps the virus latch on more tightly to human cells
- T20N May help the virus latch on tighter



Delta (B.1.617.2) Kappa (B.1.617.1) B.1.617.3

First detected: India

Key mutations:

- E484Q* May be associated with increased transmissibility; Appears to be very similar to E484K and may be associated with immune escape
- L452R May give an advantage at spreading over other variants
- P681R May help the virus latch on tighter and may result in increased transmissibility

*Absence of this mutation in the Delta variant



Epsilon (B.1.427/B.1.429)

First detected: California (USA)

Key mutations:

- L452R Reduces antibody recognition; Associated with vaccine resistance
- S13I Associated with variant that causes increased transmissibility



Zeta (P.2)

First detected: Brazil

Key mutation:

- V1176F Potential reduction in neutralization by some antibody treatments



Eta (B.1.525) B.1.1.207

First detected: United Kingdom/Nigeria

Key mutations:

- D614G Thought to make the coronavirus more infectious
- delH69V70 Alters the shape of the spike and may help it evade some antibodies
- E484K Reduces antibody recognition; Associated with vaccine resistance
- F888L Potential reduction in neutralization by some antibody treatments
- P681H Predicted to enhance systemic infection and associated with increased transmissibility



Iota (B.1.526)

First detected: New York (USA)

Key mutation:

- S477N Associated with reduced susceptibility




B.1.1.298 (Mink Variant)

Key mutation:

- D614G Thought to make the virus more infectious

We closely track literature and trends to design assays for detecting emerging SARS-CoV-2 mutations. Find the most up-to-date list of available assays at [thermofisher.com/mutationpanel](https://www.thermofisher.com/mutationpanel)


This graphic represents the most up-to-date information as of July 26, 2021.



Catalog numbers for reaction volumes

Size	No. of 10 μ L reactions (384-well plate)	No. of 20 μ L reactions (96-well plate)	Cat. No.
Small	750	374	A49785
Medium	2,500	1,250	A49786

Product	Size	Order No.
Master Mix		
TaqPath™ 1-Step RT-qPCR Master Mix, CG	5 x 1 mL	A15299
	1 x 10 mL	A15300
Recommended Controls		
AcroMetrix™ Coronavirus 2019 RNA Control	Each	954519
GeneArt™ DNA string or plasmid controls	Varies	https://www.thermofisher.com/geneart



B.1.1.7

P.1

B.1.351

Workflow

Mutation Confirmation Workflow


PCR
detection
workflow



Want to look for or confirm mutations?
Move to confirmation protocol

Use remaining RNA from prior sample extraction

RT-PCR
confirmation
workflow




Protocol using SNP assays to confirm mutations associated with SARS-CoV-2 emerging variants

Streamlined Workflow Using Gold Standard Technology


Get robust mutation detection from a simplified workflow: **1-step RT-PCR with TaqMan SNP assays**


- Use **TaqMan SNP Assays to detect mutations in RNA** extracted from SARS-CoV-2 positive samples
- Compared to typical SNP workflow this includes a **reverse-transcription (RT) step** and is done using the **Real Time setting**
- Results are analyzed in cluster plots which clearly indicate which samples have the mutation of interest and which do not



One-Step RT-PCR Using SNP Genotyping Assays

- **Start with RNA** extracted from SARS-CoV-2 positive samples
- During reverse transcription (RT), RNA is reversed transcribed to cDNA
- After reverse transcription, the PCR amplification of the cDNA template begins.
- Reaction contains a forward primer, a reverse primer, a FAM-labeled TaqMan probe, and a VIC-labeled TaqMan probe
- The fluorescence occurs only if the target sequence is complementary to the probe and amplified during PCR
- Mutation allele detected by FAM dye & reference allele by VIC dye





B.1.1.7

P.1


B.1.351

Technical Data


Reading the Results: Cluster Plot

How to interpret a cluster plot

- Cluster plots help to enable clear discrimination between reference and mutant
- Real time traces help clearly identify clusters as well as non-amplified samples
- Samples carrying reference allele will cluster along the X-axis (VIC dye; shown here in red).
- Samples carrying mutant allele will cluster along the Y-axis (FAM dye; shown here in blue).




Robust detection of S gene mutation delH69V70



of samples: 18

- delH69V70 (also known as 69-70del) is a known mutation of B.1.1.7 (501Y.V1 or UK variant)
- delH69V70 is not exclusive to B.1.1.7
- **Red**: Wild type detected
- **Blue**: delH68V70 mutation detected

Robust detection of S gene mutation N501Y



- N501Y is commonly found in:
 - B.1.1.7 (UK variant)
 - B.1.351 (South African variant)
 - P.1 (Brazilian variant)
- **Red**: Wild type detected
- **Blue**: N501Y mutation detected

of samples: 18

Thank you

Learn more and equip your lab for what's next

thermofisher.com/mutationpanel

For Research Use Only. Not for use in diagnostic procedures. © 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license.

