



Closing the Gap: Improving Diagnostics for Drug-Resistant TB


Name: Dipti Lallubhai

Date: 4 July 2024

Current state



COVID-19 changed the world. The tuberculosis (TB) response must adapt to a new reality²



A mWRD* was used as the initial diagnostic test for only 47% of the 7.5 million people newly diagnosed with TB in 2022, up from 38% in 2021 and 33% in 2020.¹



Decline in global funding available on essential TB services from U.S.\$ 6.5 billion in 2019 to U.S.\$ 5.8 billion in 2022.¹



TB deaths reduced for the first time in over a decade: **1.3 million people** died from TB in **2022**.¹

* mWRD - WHO recommended Molecular Diagnostic

1. Global tuberculosis report 2023. Geneva: World Health Organization (WHO); 2023. Accessed June 2024 [Global tuberculosis report 2023 \(who.int\)](https://www.who.int/publications/m/item/global-tuberculosis-report-2023)

2 © 2024 Cepheid. CONFIDENTIAL – SOLELY FOR INTERNAL USE BY CEPHEID ASSOCIATES.



WHO Global Report 2023

Estimated number of people who developed MDR/RR-TB (incident cases) in 2022, for countries with at least 1000 incident cases²



Three countries accounted for 42% of the estimated global number of people who developed MDR/RR-TB in 2022:

- India (27%),
- Philippines (7.5%)
- Russian Federation (7.5%)

² The eight countries ranked in descending order of the total number of RR-TB incident cases in 2022 are India, the Philippines, the Russian Federation, Indonesia, China, Pakistan, Myanmar and Nigeria.

Access to TB Diagnostics


1 62% of people have no access to rapid molecular tests for TB

2 More than **4 Mn** people are undiagnosed, untreated, and potentially transmit the disease in the community

3 More than **400k** people are infected with drug-resistant TB (MDR, RR-TB).
63% treatment success rate for people who started on 2nd line treatment

World Health Organization (WHO) 2024

Consolidated Guidelines and Operational Handbook^{*^}



- In adults with bacteriologically confirmed pulmonary TB, Xpert MTB/XDR[^] may be used on sputum for the **initial detection of resistance to isoniazid and fluoroquinolones rather than culture-based phenotypic DST.**
- In adults with bacteriologically confirmed pulmonary TB and resistance to rifampicin, Xpert MTB/XDR may be used on sputum for the initial detection of resistance to ethionamide **rather than DNA sequencing of the *inhA* promoter.**
- In adults with bacteriologically confirmed pulmonary TB and resistance to rifampicin, Xpert MTB/XDR may be used on sputum for the **initial detection of resistance to amikacin, rather than culture-based phenotypic DST.**

 **WHO recommends Xpert MTB/XDR as a low-complexity automated NAAT for the detection of resistance to isoniazid and second-line anti-TB drugs**


^{*}WHO operational handbook on tuberculosis. Module 3: diagnosis – rapid diagnostics for tuberculosis detection, third edition. Geneva: World Health Organization; 2024. Licence: CC BY-NC-SA 3.0 IGO.

[^]The WHO Operational Handbook contains information about unapproved uses of Xpert[®] MTB/XDR. Consult the instructions for use for the approved intended use.

World Health Organization (WHO)

Operational Handbook on Tuberculosis

Module 3: Diagnosis - Rapid diagnostics for tuberculosis detection, 2021 update.*



Classes of technology and associated products evaluated ^{^1}

Technology Class

Low complexity automated

Nucleic Acid Amplification Test (NAATs) for detection and resistance to isoniazid and second-line anti-TB agents

- Xpert[®] MTB/XDR^{*^} (Cepheid)

Moderate complexity automated NAATs for detection of resistance to rifampicin and isoniazid

- 7 products recommended

High complexity hybridization-based NAATs for detection of resistance to pyrazinamide

- 1 product recommended

Source

*WHO operational handbook on tuberculosis. Module 3: diagnosis – rapid diagnostics for tuberculosis detection, third edition. Geneva: World Health Organization; 2024. Licence: CC BY-NC-SA 3.0 IGO.

[^]The WHO Operational Handbook contains information about unapproved uses of Xpert[®] MTB/XDR. Consult the instructions for use for the approved intended use.

*CE-IVD. *In Vitro* Diagnostic Medical Device. May not be available in all countries. Not available in the United States.

The Cepheid Solution

The Continued Evolution of TB Diagnostics

Xpert® MTB/RIF*

The original test



Results in **110 minutes**¹

Complex technology made easy

Xpert MTB/RIF Ultra[^]

The better test



Results in **<80 minutes**²

Exceptional clinical performance

Xpert MTB/XDR[^]

Product to aid diagnosis of multidrug-resistant MTB



Results in **<90 minutes**³

Expanded diagnostic scope

1. Xpert MTB/RIF Package Insert, 301-0191,
2. Xpert MTB/RIF Ultra Package Insert, 301-5987
3. Xpert MTB-XDR ENG PI 302-3514

* CE-IVD. *In Vitro* Diagnostic Medical Device. May not be available in all countries. [^] CE-IVD. *In Vitro* Diagnostic Medical Device. May not be available in all countries. Not available in the United States.

Intended Use of the Xpert® MTB/XDR*

- The Xpert MTB/XDR test, performed on the GeneXpert Instrument Systems, is a qualitative, nested real-time polymerase chain reaction (PCR) in vitro diagnostic test for the detection of extensively drug resistant (XDR) Mycobacterium tuberculosis (MTB) complex. DNA in unprocessed sputum samples, concentrated sediments prepared from
- In specimens where MTB is detected, the Xpert MTB/XDR test can also detect **isoniazid** (INH) resistance associated mutations in the ***katG*** and ***fabG1*** genes, ***oxyR-ahpC*** **intergenic region** and ***inhA*** promoter; **ethionamide** (ETH) resistance associated with ***inhA*** promoter mutations only; **fluoroquinolone** (FLQ) resistance associated mutations in the ***gyrA*** and ***gyrB*** quinolone resistance determining regions (QRDR); and **second line injectable drug** (SLID) associated mutations in the ***rrs*** gene and the ***eis*** promoter region.
- The Xpert MTB/XDR test is intended for use as a **reflex test for a specimen (unprocessed sputum, concentrated sputum sediments, or MGIT culture) that is determined to be MTB positive.** This test is intended as an aid in the diagnosis of XDR tuberculosis (TB) when used in conjunction with clinical and other laboratory findings.




Source

1. Xpert MTB-XDR Package Insert, 302-3514
2. * CE-IVD. In Vitro Diagnostic Medical Device. May not be available in all countries. Not available in the United States

Xpert[®] MTB/XDR*

Diagnostic Pathway for Accurate Results



Fast molecular DST allows more people to start appropriate treatment on the same day

Source

1. Xpert MTB-XDR Package Insert, 302-3514
2. Xpert MTB/RIF Ultra Package Insert, 301-5987

*CE-IVD. *In Vitro* Diagnostic Medical Device. May not be available in all countries. Not available in the United States.

^US-IVD and CE-IVD. *In Vitro* Diagnostic Medical Device. May not be available in all countries

Xpert[®] MTB/XDR*

Clinical Performance

Vs. phenotypic DST

Vs. sequencing

	Vs. phenotypic DST		Vs. sequencing	
	Sensitivity	Specificity	Sensitivity	Specificity
isoniazid	91.4%	99.1%	98.8%	98.7%
fluoroquinolones	93.1%	98.5%	93.3%	100%
amikacin	91.9%	99.4%	96.4%	100%
kanamycin	87.9%	99.6%	96.7%	100%
capreomycin	84.0%	100%	96.3%	100%
ethionamide	64.7*	98.3%	97.2%	100%

➔ Sequencing as reference standard used the same gene targets

- For ethionamide, Xpert MTB/XDR* and sequencing only target the inhA promoter region, therefore have a higher discrepancy vs. phenotypic DST

Source: Xpert MTB-XDR 302-3514

* CE-IVD. In Vitro Diagnostic Medical Device. May not be available in all countries. Not available in the United States.

Xpert[®] MTB/XDR*

Clinical Performance (Prospective Specimens)

	Vs. phenotypic DST		Vs. sequencing	
	Sensitivity	Specificity	Sensitivity	Specificity
isoniazid	95.0%	95.5%	96.0%	97.7%
fluoroquinolones	94.0%	94.6%	97.1%	99.0%
amikacin	85.7%	98.4%	73.5%	99.3%
kanamycin	91.7%	92.1%	89.5%	98.4%
capreomycin	74.6%	99.4%	66.2%	99.8%
ethionamide	53.3% [^]	95.2%	96.4%	98.9%

 Clinical performance in sputum samples is comparable across different geographies in retrospective and prospective data sets.

Source

1. Xpert MTB-XDR Package Insert, 302-3514

* CE-IVD. In Vitro Diagnostic Medical Device. May not be available in all countries. Not available in the United States.

[^] Several specimens with A90V/S91P/D94A mutations in the *gyrA* gene were detected as susceptible by pDST and resistant by the test, resulting in lower specificity. Several specimens with *eis* promoter mutations and *rrs* wild type gene were detected as susceptible by pDST and resistant by the test, resulting in lower specificity. Reporting of ETH resistance is based only on the detection of *inhA* promoter mutations, resulting in a lower sensitivity

Analytical Performance of the Xpert[®] MTB/XDR* Assay

Georghio et al. Diagnostic Microbiology & Infectious Disease Journal, 2021

Highlights

➔ The Xpert MTB/XDR assay demonstrated equivalent limit of detection to Xpert MTB/RIF.

➔ The Xpert MTB/XDR assay detected 100% of tested resistance mutations and showed some utility for resistance detection in strain mixtures.

➔ For the hetero-resistance assessment, there was 100% detection of resistance when resistant populations comprised 10% of the mixture for INH resistance, 25% for FQ resistance, 50% for ETH, AMK and KAN resistance, and 60% for CAP resistance.

➔ The Xpert MTB/XDR assay reliably detects a wide range of globally relevant isoniazid, ethionamide, fluoroquinolone and second-line injectable resistance mutations.

➔ The Xpert MTB/XDR assay is a reliable, sensitive assay for tuberculosis and expanded resistance detection.

1. Georghiou S. B, et al, Analytical performance of the Xpert MTB/XDR assay for tuberculosis and expanded resistance detection, *Diagnostic Microbiology & Infectious Disease*, 2021. <https://doi.org/10.1016/j.diagmicrobio.2021.115397>

* CE-IVD. *In Vitro* Diagnostic Medical Device. May not be available in all countries. Not available in the United States.



Q&A



Thank You

www.cepheid.com