

# **STANDARD Q COVID-19 Ag Test**



---

## **Product Introduction**

**Cédric Jo / International Project Coordinator**

23 October 2020

## Contents

1. About SD BIOSENSOR
2. STANDARD Q COVID-19 Ag : Specifications
3. Control and Performance
4. Trouble shooting
5. FAQ

SD BIOSENSOR

# 01

## About SD BIOSENSOR

---



# About SD BIOSENSOR

*SD BIOSENSOR has been developing and manufacturing innovative diagnostic solutions focused on*

- **Immunoassay**

Immunoassay reagent manufacturing

- **Molecular Diagnostic Development**

Nucleic Acid Amplification reagent, POC molecular cartridge manufacturing, Nucleic Acid Extraction

- **Instruments Development**

Development of POCT system, Retention of LIS/HIS-applied technology, Lab System



**Chronic Care**  
BGMS/LipidoCare/  
MultiCare  
**Screening Test**



**STANDARD Q**  
Rapid diagnostics test




**STANDARD F**  
Fluorescent  
immunoassay



**STANDARD E**  
ELISA



**STANDARD M**  
Molecular POC system



**Confirmatory Test**

# About SD BIOSENSOR

We devote ourselves to improve human health by developing innovative products.

## 2010

**2010.12**

- Founding of SD BIOSENSOR

## 2014~2019

**2019.08 The Global Fund ERPD Approved**

- HIV/Syphilis Combo

**2019.04 UNICEF long term supply agreement signed**

- Arbo Panel I (Zika, Dengue, Chikungunya, Yellow fever)

**2016.09**

- Zika IgG/IgM

**2015.04 World 1st**

- MERS-CoV Antigen

**2014.12 WHO EUAL**

- Ebola Zaire Antigen

## 2020~

**FDA EUA in progress**

- Q COVID-19 Ag rapid
- Q COVID-19 IgM/IgG Combo rapid
- Q COVID-19 IgM/IgG Plus rapid
- F COVID-19 Ag FIA

**2020.06**

- E COVID-19 Total Ab ELISA
- Q HIV 1/2 Ab 3-Line **WHO PQ Approved**

**2020.05 WHO PQ Approved**

- Q HIV/Syphilis diagnosis kit

**2020.04 FDA EUA Approved**

- M nCoV Real-time detection kit

**2020.03**

- Q HCV Ab **WHO PQ approved**
- Q Malaria Ag **WHO PQ approved**

**CE registration**

- Q COVID-19 Ag rapid
- Q COVID-19 IgM/IgG Combo rapid
- F COVID-19 Ag FIA
- F COVID-19 IgM/IgG Combo FIA

**2020.02 CE registration /KFDA EUA Approved**

- M nCoV Real-time detection kit

# 02


## STANDARD Q COVID-19 Ag Test : Specifications

---



# ① Kit introduction

## Contents & specification



Category	Details
<b>Intended use</b>	<b>Rapid chromatographic immunoassay for qualitative detection</b> of specific SARS-CoV-2 antigen
<b>Contents (25T/kit)</b>	<ul style="list-style-type: none"><li>• Test device (individual aluminum pouch) x 25</li><li>• Sterile swab x 25</li><li>• Extraction buffer x 25</li><li>• Nozzle cap x 25</li><li>• IFU</li></ul>
<b>Sample type</b>	<b>Nasopharyngeal swab</b> <u>** Nasal swab will be added soon</u>
<b>Sample volume</b>	3 drops of mixed specimen with extraction buffer
<b>Testing time</b>	<b>15 ~ 30 minutes</b> (Do not read test results after 30 mins.)
<b>Storage temperature</b>	2~30°C (36~104°F)
<b>Operating temperature</b>	15~30°C (59~86°F) <u>** We plan to improve operating temperature until 40 °C</u>
<b>Cat. no.</b>	09COV30D (25T/kit)

## ① Kit introduction

### *Clinical Evaluation*

	Brazil	Germany	Overall
<b>Sensitivity</b> <b>(Ct ≤ 25)</b>	95.92% (47/49, 95% CI 86.02-99.50%)	100% (21/21, 95% CI 83.89-100%)	97.14% (68/70, 95% CI 90.06-99.65%)
<b>Sensitivity</b> <b>(Ct ≤ 33)</b>	91.92% (91/99, 95% CI 84.70-96.45%)	87.80% (36/41, 95% CI, 73.80-95.92%)	90.71% (127/140, 95% CI 84.64-94.96%)
<b>Sensitivity</b> <b>(0 ≤ from the symptom onset days ≤ 3)</b>	95% (19/20, 95% CI 75.13-99.87%)	85.71% (18/21, 95% CI, 63.66-96.95%)	90.24% (37/41, 95% CI 76.87 – 97.28%)
<b>Sensitivity</b> <b>(from the symptom onset days ≤ 7)</b>	90.7% (88/97, 95% CI 83.12-95.67%)	80% (28/35, 95% CI 63.06-91.56%)	87.88% (116/132, 95% CI 81.06-92.91%)
<b>Clinical Sensitivity</b>	88.68% (94/106, 95% CI 81.06-94.01%)	76.6% (36/47, 95% CI 61.97-87.70%)	84.97% (130/153, 95% CI 78.3-90.23%)
<b>Clinical Specificity</b>	97.6% (287/294, 95% CI 95.2-98.8%)	99.3% (1203/1212, 95% CI 98.6-99.6%)	98.94% (1490/1506, 95% CI 98.28-99.39%)

[https://www.finddx.org/wp-content/uploads/2020/10/SDQ-Ag-Public-Report\\_20201016-v1-1.pdf](https://www.finddx.org/wp-content/uploads/2020/10/SDQ-Ag-Public-Report_20201016-v1-1.pdf)

## ① Kit introduction

---

### *Analytical Performance*



#### **Limit of Detection (LoD)**

- The SARS-CoV-2 positive specimen was prepared by spiking Inactivated SARS-CoV-2 (2019-nCoV) NCCP 43326/2020/Korea strain to SARS-CoV-2 negative nasopharyngeal swab confirmed with PCR.
- LOD is determined as  $3.12 \times 10^{2.2}$  TCID<sub>50</sub>/ml ( $1.12 \times 10^2$  PFU/ml) for direct Nasopharyngeal swab,  $5 \times 10^{3.2}$  TCID<sub>50</sub>/ml for Nasopharyngeal swab stored in VTM by testing serially diluted the mock positive specimen.

## ② Biosafety requirements

*Material Required (Not provided)*

All the Personal Protective Equipment is disposable




Disposable container



## ② Biosafety requirements


*Recommendation of  
Centers for Disease Control and Prevention*

### [ A Biosafety Level 2 (BSL-2) Facility ]



Agents	Risk Group 2
Practices	BSL-1 plus: -Limited access -Biohazard signage -Sharps precautions -Biosafety manual
Safety Equipment	-Use of BSCs for aerosol protection -PPE-lab coats, gloves, face/eye protection
Facilities	BSL-1 plus: -Autoclave available -Directional air

### [ A Biosafety Level 3 (BSL-3) Facility ]



Agents	Risk Group 3
Practices	BSL-2 plus: -Controlled access -Decon of all waste and linens -Medical Surveillance
Safety Equipment	-Use of BSC's for all work -PPE-protective clothing, gloves, respiratory protection if needed
Facilities	BSL-2 plus: -Physical separation -Self-closing, double-door access -Negative airflow

***“STANDARD Q COVID-19 Ag Test”  
Can be used in the field***

### ※ SARS-CoV-2 in Extraction Buffer Inactivation Test

Extraction buffer	1) Virus spiking	Result	
STANDARD Q COVID-19 Extraction Buffer	O	1 minute incubation : 2)CPE	Virus Activated
	O	2 ~ 40 minutes incubation : No CPE	Virus Inactivated
	X	No CPE	Negative Control

1) SARS-CoV-2 titer :  $2.5 \times 10^{4.3}$ TCID<sub>50</sub>/mL 2) CPE : Cytopathic effect

**“ The SARS-CoV-2 virus will be inactivated “  
by extraction buffer within 2 minutes.**




## ② Biosafety requirements

### *Appropriate disinfection requirements*


✘ It is recommended to handle safely according to the recommendations of each region

Discard used sterile swab, extraction tube, test device, gloves, protective glasses, and protective suit into the disposal container



### ③ Test Handling

#### Test Kit Components (25 Tests)




※ Positive & Negative controls are provided separately (\*subject to be included in the kit soon)

### ③ Test Handling

---

#### *Storage & Transportation Requirements*

Test device & Buffer tube must be stored at ...




- Store the kit at 2-30°C / 36-86°F \*\* out of direct sunlight.
- Kit materials are stable until the expiration date printed on the outer box.
- Do not freeze the kit.

### ③ Test Kit Handling

---

#### *Sample Requirement*

Nasopharyngeal swab collected by sterile swab must be tested immediately



### ③ Test Kit Handling

---

#### *Precautions on testing assay handling*

1. Bring the kit contents and the specimens to room temperature before testing.
2. Do not re-use the test kit.
3. Do not use the test kit if the pouch is damaged or the seal is broken.
4. Do not use the extraction buffer tube of another lot.
5. Do not smoke, drink or eat while handling specimen.
6. Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.
7. Clean up spills thoroughly using an appropriate disinfectant.
8. Handle all specimens as if they contain infectious agents.
9. Observe established precautions against microbiological hazards throughout testing procedures.
10. Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and biohazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
11. Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.

## ④ Test Procedure



### Preparation

#### 1. Test device & Extraction buffer tube must check the expiry date

The front



The back



## ④ Test Procedure


### *Preparation*

#### 2. Check the test device and the desiccant


<Foil pouch>



<Test device>



<Desiccant>



## ④ Test Procedure

### *Collection of specimen*

#### 3. Collect nasopharyngeal specimen using provided sterile swab


3-1



Tilt the head back



3-2



- Tilt the patient's head back slightly and support it with your dominant hand


- Insert the swab into nasopharyngeal cavity

## ④ Test Procedure

### *Collection of specimen*


#### 3. Collect nasopharyngeal specimen using provided sterile swab

3-3



- Once swab is in location, rotate the swab
- Rotate the swab more than 5 times to saturate the swab tip

3-4



- Remove the aluminum cover of extraction buffer tube
- Remove the swab from the nasal cavity
- Insert the swab into an extraction buffer tube

## ④ Test Procedure

### *Collection of specimen*

#### 4. Extract the specimen

Mix the nasopharyngeal specimen with an extraction buffer tube

4-1




- Stir the swab more than 5 times.
- While stirring the swab, squeeze the buffer tube to extract specimen completely

4-2



- Remove the swab

4-3




- Assemble the nozzle cap tightly with the tube.

## ④ Test Procedure

### *Collection of specimen*

#### [ Specimens in transport media ]

1)



2)



#### ※ Available transport medium

Virus Transport Medium(VTM)	Recommended Storage Condition	
	2°C to 8°C	25°C
Copan UTM™ Universal Transport Media	12 hours	8 hours
BD™ Universal Viral Transport	12 hours	8 hours
STANDARD™ Transport Medium	12 hours	8 hours

- Using a micropipette, collect the 350µl of specimen from the collection cup or VTM.
- Mix the specimen with an extraction buffer.

- Assemble the nozzle cap tightly with the tube.

## ④ Test Procedure

---

### *Collection of specimen*

#### **Precautions**



If the specimen storage condition is out of instructions as below, do not use.

1. The Nasopharyngeal swab is stored in extraction buffer for more than 4 hours at  $5\pm 3^{\circ}\text{C}$  or 1 hour at  $20\pm 5^{\circ}\text{C}$ .
2. Freezing and thawing of Nasopharyngeal swab or the specimen in UTM is usable for less than 3 cycles.
3. The Nasopharyngeal swab is stored in UTM for more than 12 hours at  $5\pm 3^{\circ}\text{C}$  or 8 hours at  $20\pm 5^{\circ}\text{C}$ .

## ④ Test Procedure

---

### *Collection of specimen*

#### Criteria for rejection of specimens




- Unaccepted specimen type. Only nasopharyngeal
- Not refrigerated or frozen properly.
- Insufficient specimen volume. (Recommendation “3 drop (90 ul) ~ 4 drop (120 ul)”)
- Failure to follow specific shipping and packaging requirements.

## ④ Test Procedure

### 6. Apply the extracted specimen

6-1



- Apply 3 drops of extracted specimen to the specimen well of the test device



- Place the test device on a flat surface.
- Dispense the specimen at 90 degree angle to allow for free falling drops and avoid bubbles.

6-2



- Read the test result in 15-30 minutes


## ④ Test Procedure

### 7. Running of the assay

Read the test result in 15 - 30 minutes

#### INTERPRETATION OF TEST RESULT

\* "C" Control Line "T" Test line



Do not read test result after 30 minutes. It may show false result.  
In case of invalid, we recommend re-test

## ④ Test Procedure

---

### **Factors that affects FALSE or Invalid results**

- Concentration of specimen
- Insufficient drop of mixed buffer
- The factors that affected cross-reactivity & Interference

# 03

## QC and QA, Performance


### *STANDARD Q COVID-19 Ag Test*



# Control and Performance

---

## 1. Internal controls



- We have internal/procedural control on test device.
- The test result is valid if control line appears

## 2. Control solution

- A separate kit is available for sale (10 positive and 10 negative tablets).
- This control tablet is to be used on test device.
- Control swab will be included in the test kit in near future.

## Control and Performance

### 3. Use of software to manage

***- STANDARD PASS Mobile application can be used***

<b>Compatible product</b>	<b><u>STANDARD Q/F COVID-19 Ag Test</u></b> (It will be updated for other COVID-19 product lines soon)
<b>Target</b>	Patients who diagnosed with STANDARD Q/F COVID-19 Ag Test product (In case of negative result, Identity assurance by issuing a STANDARD Pass)
<b>Development schedule</b>	<b><u>Until 2020/12 (launch goal)</u></b>

### 4. Status of Certification & Registration of STANDARD Q COVID-19 Ag Test

<b>Product Certification</b>	<b><u>WHO EUL, CE, TGA</u></b>
<b>Countries of Registration</b>	 Morocco Tunisia Guinea Nigeria Cameroon Uganda Mozambique Botswana Republic of South Africa Angola <u>and other 8 countries</u>

**- SD BIOSENSOR distributors can support local training and A/S**

# 05

## Trouble shooting

### *STANDARD Q COVID-19 Ag Test*



## Trouble shooting

---

- Verify labeling, IFU and procedures
- Have the same operator to re-test the specimen
- Repeat blind test by another operator
- Confirm against the reference test (WHO EUL approved PCR)

# 06

## How to order

### *STANDARD Q COVID-19 Ag Test*



## How to request purchase

- **SD BIOSENSOR Homepage ([www.sdbiosensor.com/xe/covid](http://www.sdbiosensor.com/xe/covid))**



*Submit "COVID-19 Order Information"*

Title / Name / E-mail / Account / Consignee

Quantity / Incoterms / Message

- **Distributor**

- You can contact SD Biosensor distributor in your country

- **International organization**

- This kit is eligible for procurement in different platforms (ie. WHO, Global Fund, UNICEF, AMSP Etc.)

# Thank You

Cédric Jo, International Project Coordinator  
[cedric@sdbiosensor.com](mailto:cedric@sdbiosensor.com)

<http://www.sdbiosensor.com/x/covid>

**Head office** : C-4th&5th, 16, Deogyong-daero 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA

**Manufacturing site** : 74, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA

**Contact** : [sales@sdbiosensor.com](mailto:sales@sdbiosensor.com) | [www.sdbiosensor.com](http://www.sdbiosensor.com)

06  
FAQ

*STANDARD Q COVID-19 Ag Test*



## FAQs

### Q1. Should I use only nasopharyngeal swab specimen for the STANDARD Q COVID-19 Ag Test?

A. Nasopharyngeal swab and VTM(UTM) containing the specimen can be used as the sample.

Available Transport Medium

Virus Transport Medium(VTM)	Recommended Storage Condition	
	2°C to 8°C	25°C
Copan UTM™ Universal Transport Media	12 hours	8 hours
BD™ Universal Viral Transport	12 hours	8 hours
STANDARD™ Transport Medium	12 hours	8 hours

Test procedure for using VTM sample

- 1) Collect the 350ul of specimen from the VTM. Mix the specimen with an extraction buffer
- 2) Press the nozzle cap tightly onto the tube.
- 3) Apply 3 drops of extracted specimen to the specimen well of the test device.
- 4) Read the test result in 15-30 minutes

### Q2. How can I transport the specimen without VTM?

A. The specimen in an extraction buffer can be stored up to 1 hour at  $20 \pm 5^{\circ}\text{C}$  and up to 4 hours at  $5 \pm 3^{\circ}\text{C}$ .

If the specimen is not tested immediately, it is better to store the specimen at the  $-20^{\circ}\text{C}$  up to 1 cycle.  
(Do not freeze-thawing repeat)

**Q3. It is represented that the results should be interpreted between 15 – 30 mins.**

**Is there an optimum time?**

A. Interpreting the test result is available from 15 minutes.

However, it is better to read at 30 minutes as the color scale becomes more visible to interpret.

**Q4. Is there any recommendation for handling the sample and the extraction buffer?  
(E.g. inside a BSL-2-cabinet)**

A. We recommend that it is essential to follow proper infection control measures when specimen are collected from patients with a suspected CORONAVIRUS infection.

The examiner should wear an N95 respirator mask, gown, protective glasses and gloves.

Viral testing of specimens can be handled in a BSL-2 laboratory.

It is recommend that when you mixed the specimen with an extraction buffer, you have to perform the test within an hour.

**Q5. Does the extraction buffer contains any kind of substance that inactivate viable viral particles?  
I'm asking because if so, according to biosafety regulations, the test can be done outside  
a BSL-2 cabinet.**

A. The SARS-CoV-2 virus will be inactivated by extraction buffer within 2 minutes.

**※ SARS-CoV-2 in Extraction Buffer Inactivation Test**

Extraction buffer	1) Virus spiking	Result	
STANDARD Q COVID-19 Extraction Buffer	O	1 minute incubation : 2)CPE	Virus Activated
	O	2 ~ 40 minutes incubation : No CPE	Virus inactivated
	X	No CPE	Negative control

1) SARS-CoV-2 titer :  $2.5 \times 10^{4.3}$ TCID<sub>50</sub>/mL

2) CPE : Cytopathic effect

**- It's enough to inactivate the virus at least 2 minutes incubation time.**

## FAQs

---

**Q6. The transport in the buffer is a little tricky because of the nozzle cap.**

**Do they keep the swab inside of the extraction buffer (cut the ends) and then extract just before doing the test or do they transport the sample already extracted with the nozzle cap?**

A. Please perform the test at the point of care.

### Q7. Have you performed the test from Amies solution?

A. Interference has been reported in Amies solution produced by several manufacturers.

It is understood as a phenomenon that appears because the composition is slightly different for each manufacturer.