


FIND





ASLM

TB LAM technology
development pipeline

◆ Morten Ruhwald, MD, PhD
Director of TB, FIND



CURRENT AND 2ND GENERATION LAM TEST ALERE LAM, AND Fujifilm SILVAMP TB LAM





PLHIV (5 cohorts, n=1595)

	Sensitivity [95% CI]	Specificity [95% CI]
AlereLAM	34.9 [19.5 – 50.9]	95.3 [92.2 – 97.7]
FujiLam	70.7 [59.0 – 80.8]	90.9 [87.2 – 93.7]

Source: Broger et al. J Clin Invest 2020; Broger et al, PLOS Med 2020

LARGE PROSPECTIVE TRIAL DEMONSTRATE LOT TO LOT VARIABILITY




- WHO policy review cancelled
- Manufacturing issue has been resolved, test re-designed
- Relaunch expected in 2023


CHALLENGES FOR DEVELOPING LF-LAM TEST

- LAM Structure
- Sample Matrix
- LAM concentration
- Assay Design

LAM Structure



LAM Fragmentation in urine




LAM structural Integrity and stability varies greatly across samples matrix

- Lack of pre-analytical improved reagents
- Antibodies (full coverage)
 - Antigen (Cultured LAM varies from Urinary LAM)

Sample Matrixes	 Urine	 Blood	 Stool	 CSF	 Tongue Swab	 Exhale Breath Condensate	 Sputum
LAM Concentration	1 – 2000 pg/ml	1 – 2000 pg/ml	Unknown	Unknown	Unknown	100 pg-100 µg/ml	100 pg-100 µg/ml
LAM Structure	<ul style="list-style-type: none"> Delipidated Fragmentation Stability unknown 	<ul style="list-style-type: none"> Unknow Fragmentation? 	<ul style="list-style-type: none"> Less explored unknown 	<ul style="list-style-type: none"> Less explored unknown 	<ul style="list-style-type: none"> Less explored unknown 	<ul style="list-style-type: none"> Different man cap structure Truncated 	<ul style="list-style-type: none"> Intact LAM
Advantages	<ul style="list-style-type: none"> Easy to collect Potential in EPTB and pediatric TB 	<ul style="list-style-type: none"> Homogenous Potential in EPTB and pediatric TB Treatment monitoring 	<ul style="list-style-type: none"> Less explored unknown 	<ul style="list-style-type: none"> Potential in EPTB mainly TB Meningitis 	<ul style="list-style-type: none"> Easy to collect sample Treatment monitoring 	<ul style="list-style-type: none"> Easy to collect sample High LAM 	<ul style="list-style-type: none"> High LAM Treatment monitoring
Disadvantages	<ul style="list-style-type: none"> Matrix Effect Structure stability and integrity varies greatly 	<ul style="list-style-type: none"> Form Complex with HDL 	<ul style="list-style-type: none"> Less explored Matrix effects 	<ul style="list-style-type: none"> Concentration unknown 	<ul style="list-style-type: none"> Only applicable for Pulmonary TB 	<ul style="list-style-type: none"> Only applicable for Pulmonary TB 	<ul style="list-style-type: none"> Sample heterogeneity Only applicable for Pulmonary TB

3RD GEN LAM ASSAYS




THE LAM ASSAY PIPELINE

Alere LAM


- Narrow use case but saves lives



Fujifilm



SD Biosensors



- Fluorescence based test
- Under development

Boditech



- Fluorescence based test
- Under development

Abbott



- LAM concentration device
- Under development
- Visual readout

Salus



- Integrated sample preparation device
- Under development
- Visual readout

ACKNOWLEDGEMENTS



Adam Penn-Nicholson	Natalie Luong Ba	Linh Doan	Hai Yen Le	Victor Anyebe
Linh Nguyen	Lauri	Andrea Cavallini	Hoang Lyna	Vinzeigh Leukes
RUHWARD_Morten_G2022-215	Wachira Waruhu	Sandra Kik, FIND	Laure Martine	Harisha
#WorldTbDay		Dennis Walusimbi	Alex Ogwal	#YesWeCanEndTB

Harisha Ramachandraiah	Mary Gaichiri	Pamela Nabeta
Sergio Carmona	Dennis Walusimbi	Agnes Malobela
Morten Ruhwald	Alex Ogwal	Luz Villa
Adam Penn-Nicholson	Sophia Georghiou	Sunita Singh
Kavi Velen	Sandra Kik	Kavi Ramjeet
Margaretha De Vos	Victor Anyebe	Sara Frojdo
Vinzeigh Leukes	Anita Suresh	Marta Fernandez Suarez
Mikashmi Kohli	Swapna Uplekar	Aurelien Mace
Rita Szekely	Andres de la Rossa	Berra Erkosar
Andrea Cavallini	Lauri Koivula	Frida Wambui
Laure Martine	Tim Rodwell	Karishma Saran
	Rebecca Colman	Sarah-Jane Loveday
	Sacha Laurent	

Thanks to manufacturers sharing information on prototypes and strategies
 We also thank the many partners and other donors who make the work of FIND possible.

Contact: morten.ruhwald@finddx.org